

Giemsa Stain for Smears/Films - Technical Memo

SOLUTION:

Giemsa Stock Stain

500 ml

Part 1120A

1 Liter

Part 1120B

Additionally Needed:

Alcohol, Methanol Anhydrous, ACS

Part 12236

Phosphate Buffer, pH 7.0

Part 1331

For storage requirements and expiration date refer to individual bottle labels.

APPLICATION:

Newcomer Supply Giemsa Stain is a simple one-step method designed to demonstrate differential staining of cells types in peripheral blood smears and bone marrow smears/films as well as a method for detecting rickettsia, bacteria and parasites.

METHOD:

Solutions: All solutions are manufactured by Newcomer Supply, Inc.

All Newcomer Supply stain procedures are designed to be used with Coplin jars filled to 40 ml following the provided staining procedure.

STAINING PROCEDURE:

Prepare a well-made blood smear or bone marrow smear/film with a focus on uniform cell distribution.

- a. *Proceed with either the thin or thick smear/film staining method.*

Thin Smear/Film Staining Method: See Procedure Notes #1 and #2.

1. Allow smear to thoroughly air-dry prior to staining.
2. Fix smear in Methanol: 1-2 minutes.
3. Air-dry slides in a vertical position.
4. Prepare fresh 1:20 Working Giemsa Stain; combine and mix well.
 - a. *Giemsa Stock Stain* 2 ml
 - b. *Phosphate Buffer, pH 7.0 (1331)* 40 ml
5. Stain in Working Giemsa Stain for 20-30 minutes.
6. Rinse briefly in Phosphate Buffer, pH 7.0 or distilled water.
7. Air-dry slides in a vertical position.
8. If coverslip is preferred, coverslip air-dried slides with compatible mounting medium.

Thick Smear/Film Staining Method: See Procedure Notes #1 and #2.

1. Allow smear to thoroughly air-dry prior to staining.
2. Proceed directly to stain; do not place in fixative.
3. Prepare fresh 1:50 Working Giemsa Stain; combine and mix well.
 - a. *Giemsa Stock Stain* 1 ml
 - b. *Phosphate Buffer, pH 7.0 (1331)* 50 ml
4. Stain in Working Giemsa Stain for 50 minutes.
5. Rinse briefly in Phosphate Buffer, pH 7.0 or distilled water.
6. Air-dry slides in a vertical position.
7. If coverslip is preferred, coverslip air-dried slides with compatible mounting medium.

RESULTS:

Erythrocytes	Orange - pink to rose
Platelets	Red to purple granules with blue halo

Granulocytes

Neutrophils	Nucleus - Dark blue to violet
-------------	-------------------------------

Cytoplasm - Pink

Granules - Purple to lilac

Eosinophils	Nucleus - Blue
-------------	----------------

Granules - Orange to pink

Basophils	Nucleus - Deep blue to violet
-----------	-------------------------------

Granules - Deep blue to violet

Mononuclear Cells

Lymphocytes	Nuclei - Deep blue to violet
-------------	------------------------------

Cytoplasm - Light blue

Monocytes	Nuclei - Light blue/purple
-----------	----------------------------

Cytoplasm - Pale gray/blue

Mast cells	Nuclei - Deep blue to violet
------------	------------------------------

Granules - Deep blue-violet

Malarial parasites	Nucleus - Red chromatin dot
--------------------	-----------------------------

Cytoplasm - Blue

Rickettsia	Bluish purple
------------	---------------

Bacteria	Blue
----------	------

PROCEDURE NOTES:

1. Timings are suggested ranges. Optimal staining times will depend upon smear/film thickness and preference of stain intensity.
2. Smears/films containing primarily normal cell populations require minimum staining time.
 - a. *Immature cells may require a longer staining time.*
 - b. *Bone marrow smears/films may require a longer staining time.*

REFERENCES:

1. Bailey, W. Robert and Elvyn Scott. *Diagnostic Microbiology*. 4th ed. St Louis: C. V. Mosby Company, 1974. 394.
2. Garcia, Lynne Shore. *Diagnostic Medical Parasitology*. 5th ed. Washington DC: ASM Press, 2007. 888-889.
3. McPherson, Richard and Matthew Pincus. *Henry's Clinical Diagnosis and Management by Laboratory Methods*. 22nd ed. Philadelphia: Elsevier Saunders, 2011. 522-531.
4. Modifications developed by Newcomer Supply Laboratory.