

Phosphotungstic Acid Hematoxylin (PTAH) Stain - Technical Memo

SOLUTION:

Phosphotungstic Acid Hematoxylin (PTAH) Stain

500 ml

Part 1334A

1 Liter

Part 1334B

Additionally Needed:

Phosphotungstic Acid Hematoxylin (PTAH) Control Slides	Part 4565
Xylene, ACS	Part 1445
Alcohol, Ethyl Denatured, 100%	Part 10841
Alcohol, Ethyl Denatured, 95%	Part 10842
Zenker Fixative, Modified, Zinc Chloride	Part 1461
Acetic Acid, Glacial, ACS	Part 10010
Potassium Permanganate 0.25%, Aqueous	Part 133931
Oxalic Acid 5%, Aqueous	Part 1293
Coplin Jar, Plastic	Part 5184 (for microwave modification)

For storage requirements and expiration date refer to individual product labels.

APPLICATION:

The Newcomer Supply Phosphotungstic Acid Hematoxylin (PTAH) Stain procedure, with included microwave modification, is used for the demonstration of collagen, muscle striations and central nervous system (CNS) structures.

METHOD:

Fixation: Formalin 10%, Phosphate Buffered (1090)

Technique: Paraffin sections cut at 4 microns

Solutions: All solutions are manufactured by Newcomer Supply, Inc.

All Newcomer Supply stain procedures are designed to be used with Coplin jars filled to 40 ml following the provided staining procedure.

PRESTAINING PREPARATION:

- If necessary, heat dry tissue sections/slides in oven.
- Prepare Zenker Fixative Working Solution; combine and mix well.
Zenker Fixative, Modified, Zinc Chloride (1461) 38 ml
Acetic Acid, Glacial, ACS (10010) 2 ml

STAINING PROCEDURE:

- Deparaffinize sections thoroughly in three changes of xylene, 3 minutes each. Hydrate through two changes each of 100% and 95% ethyl alcohols, 10 dips each. Wash well with distilled water.
 - See Procedure Notes #1 and #2.
- Fix in Zenker Fixative Working Solution (Step #2) at 56°C; 3 hours.
Microwave Modification: See Procedure Note #3.
 - Place slides in a plastic Coplin jar containing Zenker Fixative Working Solution. Microwave 5 minutes at 60°C.
- Wash well in three changes of tap water; rinse in distilled water.
- Place in Potassium Permanganate 0.25%, Aqueous (133931) for 10 minutes.
- Wash in three changes of tap water; rinse in distilled water.
- Place in Oxalic Acid 5%, Aqueous (1293) for 10 minutes.
- Wash in three changes of tap water; rinse in distilled water.
- Place in PTAH Stain for 12-24 hours at room temperature, or 2 hours at 56°C.
 - See Procedure Note #4.
- Microwave Modification:**
 - Place slides in a plastic Coplin jar containing PTAH Stain. Microwave 7 minutes at 70°C.
- Dehydrate quickly in two changes each of 95% and 100% ethyl alcohol. Clear in three changes of xylene, 10 dips each; coverslip with compatible mounting medium.
 - Dehydrate quickly; alcohol may extract stain from sections.

RESULTS:

Collagen, cartilage, elastic fibers	Deep reddish brown
Muscle striations, fibrin, keratin	Dark blue
Glia fibers	Dark blue
Myelin	Lighter blue
Neurons	Salmon/pink
Nuclei	Blue

PROCEDURE NOTES:

- Drain slides after each step to prevent solution carry over.
- Do not allow sections to dry out at any point during procedure.
- The microwave procedure was tested using a laboratory-grade microwave oven. This procedure is a guideline and techniques should be developed for use in your laboratory.
- Adjust staining time according to preference of intensity. Suggested staining time at 37°C is 18 hours.
- If using a xylene substitute, follow manufacturer's recommendation for deparaffinization and clearing steps.

REFERENCES:

- Bancroft, John D. and Marilyn Gamble. *Theory and Practice of Histological Techniques*. 6th ed. Oxford: Churchill Livingstone Elsevier, 2008. 130-131.
- Carson, Freida L. and Christa Hladik Cappellano. *Histotechnology: A Self-Instructional Text*. 4th ed. Chicago: ASCP Press, 2015. 178-180, 201-202.
- Sheehan, Dezna C. and Barbara B. Hrapchak. *Theory and Practice of Histotechnology*. 2nd ed. St. Louis: Mosby, 1980. 193-194.
- Modifications developed by Newcomer Supply Laboratory.