

## Brown-Hopps Modified Gram Stain with Tartrazine - Technical Memo

### SOLUTION:

Tartrazine Stain 1.5%, Aqueous	<b>500 ml</b> Part 14015A	<b>1 Liter</b> Part 14015B
--------------------------------	------------------------------	-------------------------------

### Additionally Needed:

Gram, Multi-Tissue, Artificial Control Slides	Part 4256	or	Gram+ & Gram- Bacteria, Artificial Control Slides	Part 4255
Xylene, ACS	Part 1445			
Alcohol, Ethyl Denatured, 100%	Part 10841			
Alcohol, Ethyl Denatured, 95%	Part 10842			
Crystal Violet Stain 1%, Aqueous, Brown-Hopps	Part 1041			
Iodine, Gram, Aqueous	Part 1140			
Acetone, ACS	Part 10014			
Basic Fuchsin Stain 0.25%, Aqueous	Part 1011			
Gallego Solution	Part 1098			
Acetone-Xylene 1:1	Part 10015			

**For storage requirements and expiration date refer to individual product labels.**

### APPLICATION:

Newcomer Supply Brown-Hopps Modified Gram Stain with Tartrazine is a modification of the original Gram stain technique. Tartrazine, a synthetic water soluble, lemon yellow colored azo dye, provides a safe alternative to picric acid stains. Tartrazine Stain 1.5%, Aqueous provides yellow background staining, replacing the picric acid-acetone counterstain and the disposable issues associated with picric acid.

### METHOD:

**Fixation:** Formalin 10%, Phosphate Buffered (Part 1090)

**Technique:** Paraffin sections cut at 4 microns

**Solutions:** All solutions are manufactured by Newcomer Supply, Inc.

All Newcomer Supply stain procedures are designed to be used with Coplin jars filled to 40 ml following the provided staining procedure.

### STAINING PROCEDURE:

1. If necessary, heat dry tissue sections/slides in oven
2. Deparaffinize sections thoroughly in three changes of xylene, 3 minutes each. Hydrate through two changes each of 100% and 95% ethyl alcohols, 10 dips each. Wash well with distilled water.
  - a. See Procedure Notes #1 and #2.
3. Stain in Crystal Violet Stain 1%, Aqueous, Brown-Hopps (1041) for 2 minutes.
4. Rinse well in distilled water.
5. Mordant in Iodine, Gram, Aqueous (1140) for 5 minutes.
6. Rinse well in distilled water.
7. Blot excess water from slide; decolorize one slide at a time in Acetone (10014) until the blue stops running, 1-2 dips.
  - a. Sections should be very light gray in color.
8. Quickly rinse in running tap water.
9. Place in Basic Fuchsin Stain 0.25% Aqueous (1011) for 5 minutes.
10. Rinse well in running tap water.
11. Differentiate sections in Gallego Solution (1098) for 5 minutes.
12. Rinse in running tap water. Blot water off slide(s) but not to dryness.
  - a. Proceed with Steps #13 to #16 one slide at a time.
13. Place in Tartrazine Stain 1.5%, Aqueous for 1 minute.
14. Rinse well in running tap water.
15. Dip in Acetone, 1-2 quick dips.
16. Dip in Acetone-Xylene 1:1 (10015) for 5 dips.
17. Clear in three changes of xylene, 10 dips each; coverslip with compatible mounting medium.

### RESULTS:

Gram-positive bacteria	Blue
Gram-negative bacteria	Red
Nuclei	Red
Background tissue	Yellow

### PROCEDURE NOTES:

1. Drain slides after each step to prevent solution carry over.
2. Do not allow sections to dry out at any point during procedure.
3. If using a xylene substitute, closely follow the manufacturer's recommendations for deparaffinization and clearing steps.

### REFERENCES:

1. Brown, Robert C., and Howard C. Hopps. "Staining of Bacteria in Tissue Sections: A Reliable Gram Stain Method." *American Journal of Clinical Pathology* 60.2 (1973): 234-240.
2. Dapson, Janet Crookham, and Richard Dapson. *Hazardous Materials in the Histopathology Laboratory: Regulations, Risks, Handling, and Disposal*. 4th ed. Battle Creek, MI: Anatech, 2005. 150, 182, 266.
3. Mercado, Gene. "Modifying the Modification: How to Redden Shy Gram Negatives." *The Journal of Histotechnology* 25.2 (2002): 115-116.
4. Sheehan, Dezna C., and Barbara B. Hrapchak. *Theory and Practice of Histotechnology*. 2nd ed. St. Louis: Mosby, 1980. 235.
5. Modifications developed by Newcomer Supply Laboratory.