

## Wolbach Giemsa Stain - Technical Memo

### SOLUTIONS:

Giemsa Stock Stain, Wolbach

500 ml

Part 1121A

1 Liter

Part 1121B

### Additionally Needed:

Giemsa Control Slides	Part 4240
Xylene, ACS	Part 1445
Alcohol, Ethyl Denatured, 100%	Part 10841
Alcohol, Ethyl Denatured, 95%	Part 10842
Acid Alcohol 1%	Part 10011
Alcohol, Methanol Anhydrous, ACS	Part 12236
Rosin 10%, Alcoholic	Part 13398

*For storage requirements and expiration date refer to individual product labels.*

### APPLICATION:

Newcomer Supply Wolbach Giemsa Stain is used for differential staining of hematopoietic tissue and demonstration of bacteria and rickettsia that may be present in the sections

### METHOD:

**Fixation:** Recommended for hematopoietic tissue:

- Zenker Fixative, Modified, Zinc Chloride (Part 1461)
- B-5 Fixative Modified, Zinc Chloride (Part 1015)
- Formalin 10%, Phosphate Buffered (Part 1090)

**Technique:** Paraffin sections cut at 4 microns

**Solutions:** All solutions are manufactured by Newcomer Supply, Inc.

All Newcomer Supply stain procedures are designed to be used with Coplin jars filled to 40 ml following the provided staining procedure.

### PRESTAINING PREPARATION:

- If necessary, heat dry tissue sections/slides in oven.
- Prepare fresh Working Wolbach Giemsa Stain Solution; combine and mix well.
  - Distilled water 50 ml
  - Giemsa Stock Stain, Wolbach 2 ml
  - Methanol (12236) 2 ml

### STAINING PROCEDURE:

- Deparaffinize sections thoroughly in three changes of xylene, 3 minutes each. Hydrate through two changes each of 100% and 95% ethyl alcohols, 10 dips each. Wash well with distilled water.
  - See Procedure Notes #1 and #2.
- Treat slides with Acid Alcohol 1% (10011) Solution for 5 minutes.
  - See Procedure Note #3.
- Wash in running tap water for 5 minutes.
- Stain in fresh Working Wolbach Giemsa Stain Solution (Step #2) for 60 minutes at room temperature.
  - See Procedure Note #4.
- Differentiate each slide individually in Rosin 10%, Alcoholic until sections are purplish-pink; 5-10 dips. Check microscopically.
  - Rinse in 100% ethyl alcohol to stop differentiation.
- Dehydrate in two changes of 100% ethyl alcohol. Clear in three changes of xylene, 10 dips each; coverslip with compatible mounting medium.
  - Do not use a 95% ethyl alcohol dehydration step.

### RESULTS:

Nuclei	Blue/violet
Cytoplasm	Pink/rose to lighter blue shades
Bacteria	Blue
Rickettsia, inclusions	Reddish-purple

### PROCEDURE NOTES:

- Drain slides after each step to prevent solution carry over.
- Do not allow sections to dry out at any point during procedure.
- Acid Alcohol treatment ensures an acid pH and improved staining.
- For increased staining, stain in Working Giemsa Stain Solution at 60°C for 60 minutes. Additional differentiation may be required.
- The color range of cells may vary depending on fixative and degree of differentiation.
- If using a xylene substitute, closely follow the manufacturer's recommendations for deparaffinization and clearing steps.

### REFERENCES:

- Luna, Lee G. *Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology*. 3rd ed. New York: Blakiston Division, McGraw-Hill, 1968. 119-120.
- Sheehan, Dezna C, and Barbara B. Hrapchak. *Theory and Practice of Histotechnology*. 2nd ed. St. Louis: Mosby, 1980. 156-157.
- Wolbach, S. Burt, and John Todd. *The Etiology and Pathology of Typhus*. S.I.: Harvard University Press, 1922. 13-14.
- Modifications developed by Newcomer Supply Laboratory.