



Mission

Keeping it real: Our job is to help you get your job done! Everything we do is designed with your needs in mind. Our goal is to make sure your goals are met.

Keeping the lines of communication open: Our technical support staff is ready to help you solve problems, resolve concerns, and make sure that your procedures run as smoothly as possible.

Keeping prices competitive: We take pride in making sure our products arrive at your lab with the fairest price possible.

Keeping a sense of humor: Life is Good! Don't forget that. We won't. If you're having a bad day, don't worry, it's our joy to make sure your day is full of HISTOCALIFRAGILISTICXPIALIDOCIOUS!

We work hard to bring you the best histology supply services available. Please call 800-383-7799 or e-mail newly@newcomersupply.com anytime with any questions or comments.

PLACING AN ORDER

To place an order, contact customer service by phone, fax, mail or email. Our experienced staff is available and ready to help you find the histology laboratory products you need. To see our current product listings, please visit our website. Our hours of operation are Monday through Friday 8 a.m. to 5 p.m. CST. Newcomer Supply reserves the right to refuse any order.

Mail Newcomer Supply Inc.
Attn: Customer Service
2505 Parview Road
Middleton, WI 53562 USA

Phone 1-800-383-7799

Fax 1-608-831-0866

Email customerservice@newcomersupply.com

Website www.newcomersupply.com

PRICING

All prices are quoted in U.S. dollars. Newcomer Supply reserves the right to revise prices without notice. Orders will be invoiced at the price in effect on the day the order is processed. Quantity and standing order discounts are available upon request. Custom orders are subject to additional charges. Pricing does not include freight, handling charges, or taxes of any kind. Customers are responsible for paying any applicable taxes once invoiced by Newcomer Supply, or shall provide up to date and accurate tax exemption certificates.

SHIPPING

In the UNITED STATES, orders are processed and shipped usually within 3 working days from receipt of the order for In Stock inventory. *Made to Order* items are processed and shipped within 1-2 weeks. *Out of Stock* and *Custom* orders ship within 2-3 weeks. If there is an unusual delay, we will notify you. Orders are shipped F.O.B. Middleton, WI or other point of shipment pre-pay via common carrier. Third party freight accounts are accepted. If an "H" appears in the "Price" column in the website or catalog, a carrier mandated hazardous chemical charge may apply. Due to certain IATA restrictions, Newcomer Supply reserves the right to refuse to ship hazardous chemicals by air. Newcomer Supply does not accept hazardous chemical returns.

TERMS

Payment terms are net 30 days. Overdue amounts are subject to carrying charges of 1.5% per month. We accept payment by check, wire transfer, EFT, Visa, and MasterCard. We reserve the right to require prepayment by those customers that cannot meet our terms. Shipments will not be made to customers with a "Past Due" invoice pending. There is a \$5.00 handling fee per shipping box added to "each invoice. Credit memos will be honored for 6 months. No cash refunds.

RETURNS

Contact Newcomer Supply prior to returning any product. If you have a problem with a product or shipment, call us within 5 business days of the issue. Returns without a Return Authorization "RA" will not be accepted. Returns must be packaged as received and delivered undamaged within two weeks of RA date. If you ordered a product in error and request a return, a 25% restocking fee is applicable. Customers are responsible for return shipment costs. Original freight costs will not be credited. All custom orders are final, no returns. Newcomer Supply does not accept hazardous chemical returns.

GENERAL INFORMATION

All products represented in the Newcomer Supply catalog are intended only for in vitro laboratory use by qualified professionals. They are not intended for drug, cosmetic, food, agricultural, pesticide, or household use. Any use other than that specifically authorized voids any and all warranties.

SAFETY INFORMATION

For shelf life and storage temperature refer to individual product label. Chemicals should be handled only by qualified individuals trained in laboratory procedures and familiar with the potential hazards. Safety Data Sheets are available on our website, by email or mail. Technical Memos describing procedural methods are available for "our customers via our website, mail or email.

LIMITED WARRANTY

NEWCOMER SUPPLY, INC. ("NEWCOMER") warrants only that products manufactured by NEWCOMER are as represented in the catalog, and that products sold under the NEWCOMER label meet the minimum published specifications as to chemical and physical properties at the time of sale. NEWCOMER is not responsible if age or conditions of shipment or storage alter the physical and/or chemical properties of any products. If any product does not meet the standards of this warranty, NEWCOMER SUPPLY will, at its sole option, either (A) replace the product or products that do not meet this warranty; or (B) refund the purchase price for the product or products that do not meet this warranty; *PROVIDED* the product or products which do not meet this warranty are returned to NEWCOMER SUPPLY, INC. for inspection within thirty (30) days of the Date of Invoice, and such inspection reveals, to the sole satisfaction of NEWCOMER, that the warranty has been breached. This warranty is exclusive. NEWCOMER specifically disclaims any and all other warranties, express or implied, including the implied warranties of merchantability, fitness for a particular purpose, non-infringement and otherwise. In no event will NEWCOMER be liable for any lost revenue, lost profit, lost data, or attorney fees, or for any special, indirect, consequential, incidental, or punitive damages, however caused and regardless of the theory of liability. In no event shall NEWCOMER'S liability, if any, whether in contract, tort (including negligence), or otherwise, exceed the purchase price actually paid for the product or products that fail to meet this warranty.

TABLE OF CONTENTS

POSITIVE CONTROL SLIDES

PG 6

Custom	14
Histopathology	6
Immunopathology	10

STAINS & STAIN KITS

PG 16

Bluing Agents	17
Buffers	25
Eosins	16
Hematoxylin	16
Special Stains	17
Stain Kits	23

CHEMICALS & REAGENTS

PG 26

Alcohols	26
Chemicals	32
Decalcifying Solutions	27
Fixatives	28
Michel Transport Medium	30
Xylenes	31

SLIDES & COVER SLIPS

PG 38

Adhesive / Charged	38
Adhesive Solutions	42
Cover Slips	40
Mounting Mediums	41
Non-Adhesive	38
Red Control Box	40

LAB ESSENTIALS & EQUIPMENT

PG 44

Essentials

Brushes	44
Microtome & Cryostat Aligners	44
Microtome Blades	44
OCT	45
PAP Pens	45
Parafilm M	46
Pens, Marker	46
Pipettes	46
Slide Boxes	47
Slide Draining Trays	48
Slide Mailers	48
Slide Presentation Folders	48
Slide Staining Trays	49
Specimen Chucks	50
Staining Dishes & Racks	50
Sterile Instrument Packs	52
Storage Systems	52
Timers	52
Xylene Pads	53

Equipment

Alcohol Recycling	53
Centrifuges	54
Cryo Bath	55
Cryostat Freeze Plate	55
Keyboards & Mouse	55
Microtomes	56
Paraffin Dispenser	56
Paraffin Wax Trimmers	56
Slide Warmers	56
Stirring Hot Plate	56
Tissue Embedding Center	57
Water Baths	57

**GROSSING & EMBEDDING****PG 58**

Base Molds	58
Baskets for Tissue Cassettes	61
Biopsy Sponge Pads & Wraps	60
Cassettes	58
Dissecting Mat	61
Dissecting Needles	67
Forceps	61
Forceps Jars	63
Hammer	67
Instrument Holder	67
Paraffin & Repellent	66
Rulers for Pathology	65
Scalpels	63
Scissors	64
Telescopic Glasses	66
Tissue Dyes & Accessories	65
Tissue Tampers	58

SAFETY**PG 68**

Air Purifier	72
Disinfectants	68
Formalin Absorbing Pads	69
Formalin Cube Containment	69
Formalin Disposal	69
Gloves Nitrile Powder Free	71
Safety Glasses	71
Sharps Containers	72
Spill Kits & Components	70

IHC**PG 74**

Buffers	74
Retrieval Solutions	75

(We are a distributor of Innovex Biosciences and Golden Bridge International Labs IHC product lines; Detection Kits, Chromogens, Antibodies and Specialty Products.)

Their part numbers and pricing are available on our website: www.newcomersupply.com

POSITIVE CONTROL SLIDES

PRODUCT INFORMATION

We have one of the largest selections of both histopathology and immunopathology control slides available. Also, we can accommodate custom requests, such as tissue placement and specialty microscope slide glass. Additional charges may apply.

All control slides are on Superfrost® Plus

All controls are cut in serial sections.

10/set = 10 unstained slides and 1 stained slide.

98/set = 98 unstained slides and 2 stained slides.

Please note the following: In human tissue unless otherwise indicated. All tissues are FFPE. All slides are heat dried. Histopathology slides are cut at 5 microns, unless otherwise specified. Standard tissue placement is mid-range of slide. Histopathology control slides are not quality controlled for IHC methods, as such, we cannot guarantee antigenicity. The intended use of these control slides is to qualitatively assess the performance of the special stain. These controls are not intended for quantitative purposes. Review the supplied stained slides with your pathologist prior to using the unstained slides to ensure that the intended analyte of interest is acceptable.

Tissue Placement

Standard tissue placement on our histology and immunopathology control slides is in the middle of the slide.

Alternately, you can request HIGH or LOW tissue placement on the slide.

Heat Drying

All Newcomer Supply control slides are dried horizontally, in an incubator at 65°C for 90 minutes. Customers may also specify how slides are oriented during drying or for no drying at all.

PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
Histopathology				
Acid Fast Bacteria (AFB), Human (Inquire for availability)	(+) Soft Tissue / (-) Lung	4010A	10/set	\$45.00
		4010B	98/set	\$323.00
Acid Fast Bacteria (AFB), Animal	(+) Goat Intestine / (-) Lung	4011A	10/set	\$45.00
		4011B	98/set	\$323.00
Alcian Blue pH 2.5, Umbilical Cord	(+) Umbilical Cord	4020A	10/set	\$45.00
		4020B	98/set	\$323.00
Alcian Blue pH 2.5, Barrett's Esophagus Limit 3 sets per order.	(+) Barrett's Esophagus / (-) Liver	4023A	10/set	\$45.00
Alcian Blue pH 2.5, Multi-Tissue	(+) Umbilical Cord	4021AMT	10/set	\$63.00
	(+) Small Intestine / (-) Liver	4021BMT	98/set	\$445.00
Alcian Blue/PAS	(+) Small Intestine	4022A	10/set	\$45.00
		4022B	98/set	\$323.00
Amoeba	(+) Intestine	4025A	10/set	\$45.00
		4025B	98/set	\$323.00
Amyloid, Animal	(+) Animal Organ	4031A	10/set	\$45.00
		4031B	98/set	\$323.00
Argentaffin	(+) Small Intestine	4035A	10/set	\$45.00
		4035B	98/set	\$323.00



PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
Argyrophil	(+) Small Intestine	4040A	10/set	\$45.00
Basement Membrane	(+) Kidney	4055A	10/set	\$45.00
		4055B	98/set	\$323.00
Bielschowsky	(+) Brain	4058A	10/set	\$45.00
		4058B	98/set	\$323.00
Bile	(+) Liver	4060A	10/set	\$45.00
		4060B	98/set	\$323.00
Blastomyces sp. (Inquire for availability)	(+) Lung / (-) Lung	4090A	10/set	\$45.00
		4090B	98/set	\$323.00
Blastomyces sp., Animal	(+) Animal Lymph Node / (-) Lung	4091A	10/set	\$45.00
		4091B	98/set	\$323.00
Calcium	(+) Placenta	4100A	10/set	\$45.00
		4100B	98/set	\$323.00
Cat Scratch Bartonella sp., Artificial	(+) Rat Lung / (-) Lung	4110A	10/set	\$45.00
		4110B	98/set	\$323.00
Coccidioides sp. (Inquire for availability)	(+) Lung / (-) Lung	4125A	10/set	\$45.00
Colloidal Iron	(+) Small Intestine	4127A	10/set	\$45.00
		4127B	98/set	\$323.00
Colloidal Iron, Multi-Tissue	(+) Small Intestine (+) Iron, Animal	4128AMT	10/set	\$63.00
		4128BMT	98/set	\$445.00
Copper, Animal	(+) Animal Liver	4130A	10/set	\$45.00
		4130B	98/set	\$323.00
Cryptococcus sp.	(+) Lung / (-) Lung	4135A	10/set	\$45.00
		4135B	98/set	\$323.00
Elastic, Aorta	(+) Aorta	4194A	10/set	\$45.00
		4194B	98/set	\$323.00
Elastic, Skin	(+) Skin	4195A	10/set	\$45.00
		4195B	98/set	\$323.00
Fite, Norcardia sp., Artificial Nocardia sp. as procedural control for leprosy.	(+) Rat Lung / (-) Lung	4215A	10/set	\$45.00
		4215B	98/set	\$323.00
Fungus, GMS, Aspergillus sp., Artificial	(+) Rat Lung / (-) Lung	4220A	10/set	\$45.00
		4220B	98/set	\$323.00
Fungus, GMS, Aspergillus sp., Animal	(+) Duck Lung / (-) Lung	4221A	10/set	\$45.00
		4221B	98/set	\$323.00
Fungus, GMS, Candida sp., Artificial	(+) Rat Lung / (-) Lung	4230A	10/set	\$45.00
		4230B	98/set	\$323.00

POSITIVE CONTROL SLIDES

PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
Histopathology				
Fungus, PAS, Aspergillus sp., Artificial	(+) Rat Lung / (-) Lung	4232A	10/set	\$45.00
		4232B	98/set	\$323.00
Fungus, PAS, Candida sp., Artificial	(+) Rat Lung / (-) Lung	4233A	10/set	\$45.00
		4233B	98/set	\$323.00
Fungus, GMS, Multi-Tissue, Artificial	(+) Aspergillus Rat Lung (+) Candida Rat Lung (+) Histoplasma Animal Infection / (-) Lung	4235AMT	10/set	\$70.00
		4235BMT	98/set	\$492.00
Giemsa	(+) Spleen	4240A	10/set	\$45.00
		4240B	98/set	\$323.00
Gram, Gram Positive & Gram Negative Bacteria, Artificial	(+) Rat Lung / (-) Lung	4255A	10/set	\$45.00
		4255B	98/set	\$323.00
Gram, Multi-Tissue, Artificial	(+) Gram Pos Rat Lung (+) Gram Neg Rat Lung / (-) Lung	4256A	10/set	\$45.00
		4256B	98/set	\$323.00
Gram, Gram Positive & Gram Negative Bacteria (Inquire for availability)	(+) Soft Tissue / (-) Lung	4257A	10/set	\$45.00
Helicobacter sp., Animal	(+) Primate Stomach / (-) Lung	4274A	10/set	\$62.00
		4274B	98/set	\$427.00
Helicobacter sp., Artificial	(+) Rat Lung / (-) Lung	4275A	10/set	\$45.00
		4275B	98/set	\$323.00
Helicobacter sp.	(+) Gastrointestinal Various / (-) Lung	4276A	10/set	\$62.00
		4276B	98/set	\$427.00
Hepatitis Limit of 3 sets per order.	(+) Liver	4280A	10/set	\$45.00
Histoplasma sp., Animal	(+) Animal Organ / (-) Lung	4290A	10/set	\$45.00
		4290B	98/set	\$323.00
Iron	(+) Spleen, Organ Varies	4320A	10/set	\$45.00
		4320B	98/set	\$323.00
Iron, Animal	(+) Animal Spleen, Organ Varies	4321A	10/set	\$45.00
		4321B	98/set	\$323.00
Legionella sp.	(+) Lung / (-) Lung	4400A	10/set	\$45.00
		4400B	98/set	\$323.00
Leptospira sp., Animal	(+) Hamster Liver / (-) Lung	4402A	10/set	\$45.00
		4402B	98/set	\$323.00



PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
Luxol Fast Blue (LFB) - Cresyl Violet	(+) Spinal Cord	4407A	10/set	\$45.00
		4407B	98/set	\$323.00
Mast Cell, Animal	(+) Animal Tumor	4412A	10/set	\$45.00
		4412B	98/set	\$323.00
Melanin	(+) Skin	4430A	10/set	\$45.00
		4430B	98/set	\$323.00
Movat-Russell Pentachrome	(+) Small Intestine (+) Lung	4440A	10/set	\$45.00
		4440B	98/set	\$323.00
Mucin Mucicarmine	(+) Colon	4455A	10/set	\$45.00
		4455B	98/set	\$323.00
Periodic Acid Schiff (PAS) Glycogen With and without amalyse digestion.	(+) Liver	4540A	10/set	\$45.00
		4540B	98/set	\$323.00
Phosphotungstic Acid Hematoxylin (PTAH)	(+) Striated Muscle	4565A	10/set	\$45.00
		4565B	98/set	\$323.00
Pneumocystis Carinii (Inquire for availability)	(+) Lung / (-) Lung	4555A	10/set	\$45.00
		4555B	98/set	\$323.00
Pneumocystis sp., Animal	(+) Primate Lung / (-) Lung	4556A	10/set	\$45.00
		4556B	98/set	\$323.00
Reticulum	(+) Liver	4620A	10/set	\$45.00
		4620B	98/set	\$323.00
Spirochete, Artificial	(+) Rat Lung / (-) Lung	4655A	10/set	\$45.00
		4655B	98/set	\$323.00
Spirochete, Animal	(+) Dog Organ / (-) Lung	4660A	10/set	\$45.00
		4660B	98/set	\$323.00
Trichrome, Liver	(+) Liver	4690A	10/set	\$45.00
		4690B	98/set	\$323.00
Trichrome, Kidney	(+) Kidney	4691A	10/set	\$45.00
		4691B	98/set	\$323.00
Trichrome, Uterus	(+) Uterus	4695A	10/set	\$45.00
		4695B	98/set	\$323.00
Trichrome, Multi-Tissue	(+) Liver (+) Kidney (+) Uterus	4693AMT	10/set	\$70.00
		4693BMT	98/set	\$492.00
Urates Limit of 3 sets per order.	(+) Gouty Tophi	4700A	10/set	\$45.00

POSITIVE CONTROL SLIDES

PRODUCT INFORMATION

We have one of the largest selections of both histopathology and immunopathology control slides available. Also, we can accommodate custom requests, such as tissue placement and specialty microscope slide glass. Additional charges may apply.

All control slides are on Superfrost® Plus

All controls are cut in serial sections.

10/set = 10 unstained slides and 1 stained slide.

98/set = 98 unstained slides and 2 stained slides.

Please note the following: Allow 2-3 weeks for delivery of large sets. In human tissue unless otherwise indicated. All tissues are FFPE. All slides are heat dried. Immunopathology slides are cut at 4 microns. Standard tissue placement is mid-range of slide. Visit our website for the most current control slide offering. The intended use of these control slides is to qualitatively assess the performance of the Immunohistochemical stain. These controls are not intended for quantitative purposes. Review the supplied stained slides with your pathologist prior to using the unstained slides to ensure that the intended analyte of interest is acceptable.

Tissue Placement

Standard tissue placement on our histology and immunopathology control slides is in the middle of the slide. Alternately, you can request HIGH or LOW tissue placement on the slide.

Heat Drying

All Newcomer Supply control slides are dried horizontally, in an incubator at 65°C for 90 minutes. Customers may also specify how slides are oriented during drying or for no drying at all.

PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
Immunopathology				
Actin, Muscle Specific	(+) Uterus	3002A	10/set	\$53.00
		3002B	98/set	\$411.00
Actin, Smooth Muscle	(+) Uterus	3003A	10/set	\$53.00
		3003B	98/set	\$411.00
Alpha Fetoprotein (AFP)	(+) Fetal Liver / (-) Liver	3055A	10/set	\$53.00
		3055B	98/set	\$411.00
BCL-2 Oncoprotein	(+) Lymphoma	3052A	10/set	\$53.00
CA-125	(+) Endometrial Carcinoma /	3140A	10/set	\$53.00
	(-) Kidney	3140B	98/set	\$411.00
Calponin	(+) Uterus	3148A	10/set	\$53.00
		3148B	98/set	\$411.00
Calretinin Limit of 3 sets per order.	(+) Mesothelioma / (-) Kidney	3150A	10/set	\$53.00
Carcinoembryonic Antigen (CEA)	(+) Colon Adenocarcinoma /	3160A	10/set	\$53.00
	(-) Kidney	3160B	98/set	\$411.00
Cat Scratch Bartonella sp., Artificial	(+) Rat Lung / (-) Lung	3165A	10/set	\$53.00
		3165B	98/set	\$411.00



PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
CD3, T-Cell	(+) Tonsil	3770A	10/set	\$53.00
		3770B	98/set	\$411.00
CD15	(+) Hodgkin's Lymphoma / (-) Myometrium	3167A	10/set	\$53.00
		3167B	98/set	\$411.00
CD20, B-Cell	(+) Tonsil / (-) Kidney	3050A	10/set	\$53.00
		3050B	98/set	\$411.00
CD23	(+) Tonsil / (-) Kidney	3173A	10/set	\$53.00
		3173B	98/set	\$411.00
CD30	(+) Hodgkin's Lymphoma / (-) Kidney	3175A	10/set	\$53.00
		3175B	98/set	\$411.00
CD31, PECAM-1	(+) Angiosarcoma (+) Tonsil	3174A	10/set	\$53.00
		3174B	98/set	\$411.00
CD34	(+) Lung (Spindle Cell Sarcoma) (+) Tonsil	3176A	10/set	\$53.00
		3176B	98/set	\$411.00
CD45, Leucocyte Common Antigen (LCA)	(+) Tonsil / (-) Adipose	3545A	10/set	\$53.00
		3545B	98/set	\$411.00
CD56	(+) Pancreas / (-) Lung	3184A	10/set	\$53.00
		3184B	98/set	\$411.00
CD117, c-Kit	(+) GIST / (-) Adipose	3170A	10/set	\$53.00
		3170B	98/set	\$411.00
Chromogranin A	(+) Pancreas / (-) Kidney	3190A	10/set	\$53.00
		3190B	98/set	\$411.00
CK5/14 + p63 + P504S, Multi-Tissue (Referred to as PIN-4 Cocktail)	(+) Prostate Adenocarcinoma (+) Benign Prostate / (-) Myometrium	3650AMT	10/set	\$121.00
		3650BMT	98/set	\$662.00
c-Myc	(+) Breast Carcinoma / (-) Lung	3171A	10/set	\$53.00
		3171B	98/set	\$411.00
Cytokeratin AE1/AE3	(+) Skin (+) Colon / (-) Myometrium	3235A	10/set	\$53.00
		3235B	98/set	\$411.00
Cytokeratin 7 (CK7)	(+) Lung Carcinoma / (-) Myometrium	3232A	10/set	\$53.00
		3232B	98/set	\$411.00
Cytokeratin 20 (CK20)	(+) Colon Adenocarcinoma / (-) Myometrium	3233A	10/set	\$53.00
		3233B	98/set	\$411.00
Cytokeratin, HMW	(+) Skin / (-) Myometrium	3237A	10/set	\$53.00
		3237B	98/set	\$411.00
Cytomegalovirus (CMV)	(+) Lung / (-) Tonsil	3240A	10/set	\$53.00
		3240B	98/set	\$411.00

POSITIVE CONTROL SLIDES

PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
Immunopathology				
Desmin	(+) Leiomyoma Uterus /	3255A	10/set	\$53.00
	(-) Tonsil	3255B	98/set	\$411.00
DOG1	(+) GIST / (-) Kidney	3257A	10/set	\$53.00
		3257B	98/set	\$411.00
E-cadherin	(+) Ductal Breast Carcinoma /	3283A	10/set	\$53.00
	(-) Myometrium			
Epithelial Membrane Antigen (EMA)	(+) Breast Carcinoma /	3285A	10/set	\$53.00
	(-) Thyroid	3285B	98/set	\$411.00
Glial Fibrillary Acidic Protein (GFAP)	(+) Glioblastoma /	3320A	10/set	\$53.00
	(-) Myometrium	3320B	98/set	\$411.00
Helicobacter sp.	(+) Gastrointestinal / (-) Lung	3276A	10/set	\$73.00
		3276B	98/set	\$492.00
Helicobacter sp., Animal	(+) Primate Stomach / (-) Lung	3274A	10/set	\$73.00
		3274B	98/set	\$492.00
Helicobacter sp., Artificial	(+) Rat Lung / (-) Lung	3270A	10/set	\$53.00
		3270B	98/set	\$411.00
Hepatitis B Core Antigen (HBcAG) Limit of 3 sets per order.	(+) Liver / (-) Kidney	3360A	10/set	\$53.00
Hepatitis B Surface Antigen (HBsAG) Limit of 3 sets per order.	(+) Liver / (-) Kidney	3365A	10/set	\$53.00
Herpes Simplex Virus, Type 1&2 (HSV-1&2) (Inquire for availability)	(+) Liver / (-) Lung	3281A	10/set	\$53.00
HMB-45	(+) Melanoma / (-) Lung	3370A	10/set	\$53.00
		3370B	98/set	\$411.00
Kappa	(+) Tonsil	3480A	10/set	\$53.00
		3480B	98/set	\$411.00
Ki-67	(+) Breast Carcinoma /	3485A	10/set	\$53.00
	(-) Adipose	3485B	98/set	\$411.00
Lambda	(+) Tonsil	3530A	10/set	\$53.00
		3530B	98/set	\$411.00
Melan-A (MART-1)	(+) Melanoma / (-) Lung	3580A	10/set	\$53.00
		3580B	98/set	\$411.00
Neuron-Specific Enolase (NSE)	(+) Pancreas / (-) Adipose	3611A	10/set	\$53.00
		3611B	98/set	\$411.00



PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
p16	(+) Breast Carcinoma / (-) Adipose	3625A	10/set	\$53.00
		3625B	98/set	\$411.00
p40	(+) Skin Squamous Cell Carcinoma / (-) Kidney	3635A	10/set	\$53.00
		3635B	98/set	\$411.00
p53	(+) Breast Carcinoma / (-) Adipose	3640A	10/set	\$53.00
		3640B	98/set	\$411.00
p63	(+) Benign Prostate / (-) Myometrium	3641A	10/set	\$53.00
		3641B	98/set	\$411.00
Pneumocystis Carinii (Inquire for availability)	(+) Lung / (-) Lung	3670A	10/set	\$53.00
		3670B	98/set	\$411.00
Prostate Specific Antigen (PSA)	(+) Prostate / (-) Lung	3681A	10/set	\$53.00
		3681B	98/set	\$411.00
S-100	(+) Melanoma / (-) Lung	3740A	10/set	\$53.00
		3740B	98/set	\$411.00
Spirochete Treponema sp., Artificial	(+) Rat Lung / (-) Lung	3746A	10/set	\$53.00
		3746B	98/set	\$411.00
Spirochete, Animal	(+) Dog Organ / (-) Lung	3745A	10/set	\$53.00
		3745B	98/set	\$411.00
Synaptophysin	(+) Pancreas / (-) Kidney	3750A	10/set	\$53.00
		3750B	98/set	\$411.00
Thyroid Transcription Factor (TTF-1)	(+) Thyroid / (-) Myometrium	3780A	10/set	\$53.00
		3780B	98/set	\$411.00
Toxoplasma sp., Animal	(+) Rat Tissue, Various / (-) Myometrium	3801A	10/set	\$53.00
		3801B	98/set	\$411.00
Uro-2™ (CK20 + p53), Multi-Tissue	(+) Colon Carcinoma (+) Benign Colon / (-) Lung	3830AMT	10/set	\$118.00
		3830BMT	98/set	\$645.00
Vimentin	(+) Melanoma / (-) Adipose	3870A	10/set	\$53.00
		3870B	98/set	\$411.00
Wilms' Tumor (WT1)	(+) Kidney Tumor / (-) Lung	3900A	10/set	\$53.00
		3900B	98/set	\$411.00

POSITIVE CONTROL SLIDES

PRODUCT INFORMATION

We offer a variety of different slide types for you to choose from. If you would like something different than our standard Superfrost® Plus offering, please let us know.

Tissue Placement

Standard tissue placement on our histology and immunopathology control slides is in the middle of the slide.

Alternately, you can request HIGH or LOW tissue placement on the slide.

Heat Drying

All Newcomer Supply control slides are dried horizontally, in an incubator at 65°C for 90 minutes. Customers may also specify how slides are oriented during drying or for no drying at all.

PRODUCT NAME	TISSUE	PART	VOLUME	PRICE
Custom				
Red Control Box				
Histopathology Quality Control Slides		*4###ARB	10/Set	\$49.00
Histopathology Quality Control Slides		*4###BRB	98/Set	\$379.00
Helicobacter, Human (Histo) Control Slides	(+) Gut Various / (-) Lung	4276ARB	10/set	\$65.00
Helicobacter, Human (Histo) Control Slides		4276BRB	98/set	\$474.00
Immunopathology Quality Control Slides		*3###ARB	10/set	\$59.00
Immunopathology Quality Control Slides		*3###BRB	98/set	\$441.00
Helicobacter, Human (Immuno) Control Slides	(+) Gut Various / (-) Lung	3276ARB	10/set	\$75.00
Helicobacter, Human (Immuno) Control Slides		3276BRB	98/set	\$526.00



PRODUCT NAME

Custom Tissue Slides

Newcomer Supply Custom Tissue Slides line features an inventory of human and animal tissue slides for your custom histochemical and immunohistochemical (IHC) staining (testing) needs. Custom tissue slides are from formalin fixed paraffin embedded (FFPE) samples and are available from a wide range of organs and a variety of disease states. Our Custom Tissue Slides are sold in sets of 10 unstained slides with a complimentary hematoxylin and eosin (H&E) stained slide included in each set.

Please call us to inquire about part number and pricing for this product line.

**NEW**

STAINS & STAIN KITS

PRODUCT INFORMATION

- Shelf Life is one year from date manufactured unless otherwise stated.
- Store at 15-30°C unless otherwise stated.

H = Hazardous

PRODUCT NAME	PART	VOLUME	PRICE	
Hematoxylin & Eosin Staining				
Hematoxylin				
Hematoxylin Stain, Harris Modified Harris without the mercury. Acetic enhanced.	1201A	500 ml	\$34.00	
	1201B	1 L	\$53.50	
	1201C	1 gal	\$121.60	
Hematoxylin Stain, Harris Non-Mercuric Stronger. Crisper.	12013A	500 ml	\$36.30	
	12013B	1 L	\$57.10	
	12013C	1 gal	\$129.80	
Hematoxylin Stain, Gill				
Gill I - Use: cytology	1180A	500 ml	\$23.70	
	1180C	1 gal	\$95.00	
Gill II - Use: cytology and histology	1180D	500 ml	\$29.90	
	1180F	1 gal	\$120.70	
Gill III - Use: histology	1180G	500 ml	\$24.70	
	1180I	1 gal	\$148.80	
Hematoxylin Stain, Mayer Modified No Chloral Hydrate	1202A	500 ml	\$34.40	
	1202B	1 L	\$54.10	
	1202C	1 gal	\$151.00	
Hematoxylin, Powder	1160A	10 g	\$38.30	
	1160B	25 g	\$105.50	
	1160C	100 g	\$161.00	
Hematoxylin Stain Set, Weigert Iron (Use: Mucin and Trichromes)				
Solution A: Ferric Chloride, Aqueous				
Solution B: Hematoxylin 1%, Alcoholic				
Each bottle 250 ml	1409B	500 ml set	\$82.80	H
Each bottle 500 ml	1409A	1 L set	\$110.80	H
Eosins				
Eosin-Phloxine Stain Set Solution A: Eosin Y Stock Stain 1%, Aqueous, 1 L Solution B: Phloxine B Stock Stain 1%, Aqueous, 100 ml	1082A	set	\$52.40	
Eosin Y Stock Stain 1%, Alcoholic (Use: Counterstain for H&E Stain)	1070B	500 ml	\$16.30	H
	1070C	1 L	\$27.40	H
	1070D	1 gal	\$66.80	H



PRODUCT NAME	PART	VOLUME	PRICE	
Eosin Y Stock Stain 1%, Aqueous (Use: Counterstain for H&E Stain)	1080B	500 ml	\$12.60	
	1080C	1 L	\$17.30	
	1080D	1 gal	\$57.10	
Eosin Y Working Solution (Use: Counterstain for H&E Stain)	1072A	500 ml	\$15.70	H
	1072B	1 L	\$25.90	H
	1072C	1 gal	\$49.30	H
Eosin Y, Powder	1060A	25 g	\$16.60	
	1060B	100 g	\$53.70	

Bluing Agents

Scott Tap Water Substitute Also see: Lithium Carbonate and Ammonium Hydroxide in the Chemicals & Reagents Section	1380A	500 ml	\$14.10	
	1380B	1 L	\$23.20	
	1380C	1 gal	\$32.80	

Special Stains

Alcian Blue 8GX, Powder	1002A	5 g	\$51.60	
Alcian Blue Stain 1%, pH 1.0 Aqueous (Use: Sulfated Acidic Mucosubstances) Shelf Life: 3 mos.	1005A	250 ml	\$43.80	
	1005B	500 ml	\$93.30	
Alcian Blue Stain 1%, pH 2.5 Aqueous (Use: Acidic Mucins & Carboxylated Mucosubstances; including Barrett's Esophageal Biopsy) Shelf Life: 6 mos.	1003A	250 ml	\$68.80	
	1003B	500 ml	\$137.40	
Alcian Blue Stain, Alcoholic (Use: Sulfated Alcian Blue for Amyloid)	1004A	500 ml	\$122.50	H
Aniline Blue, Powder	10073A	10 g	\$11.40	
	10073B	25 g	\$27.10	
	10073C	50 g	\$54.30	
Aniline Blue Stain, Aqueous (Use: Masson or McLetchie Trichrome Stain)	10072B	250 ml	\$28.90	
	10072C	500 ml	\$43.50	
Auramine O, Powder (Use: Fluorescent dye)	10088A	25 g	\$32.30	
	10088B	100 g	\$85.90	
Basic Fuchsin, Powder	1010A	10 g	\$24.80	
	1010B	25 g	\$39.00	
	1010C	100 g	\$140.30	
Basic Fuchsin Stain 0.25%, Aqueous (Use: Brown-Brenn mod. Gram Stain; Can be used for Brown-Hopps)	1011A	250 ml	\$26.10	
	1011B	500 ml	\$40.10	

STAINS & STAIN KITS

PRODUCT NAME	PART	VOLUME	PRICE
Special Stains			
Basic Fuchsin Stain 0.5%, Aqueous (Use: Brown-Hopps mod. Gram Stain)	1012A	250 ml	\$26.80
	1012B	500 ml	\$39.70
Biebrich Scarlet-Acid Fuchsin Stain, Aqueous (Use: Masson or McLetchie Trichrome Stain)	10161B	250 ml	\$45.60
	10161C	500 ml	\$76.10
Biebrich Scarlet-Acid Fuchsin Stain, Elastic-Trichrome, Aqueous (Use: Modified Elastic Trichrome)	1016A	250 ml	\$57.90
	1016B	500 ml	\$81.70
Carbol Fuchsin Stain, Ziehl-Neelson (Use: AFB Stain for Ziehl-Neelson, Fite) Shelf Life: 2 yrs.	1030A	250 ml	\$36.70
	1030B	500 ml	\$56.20
	1030C	1 L	\$96.40
Carbol Fuchsin Stain, Kinyoun (Use: Kinyoun for AFB) Shelf Life: 6 mos.	1031A	250 ml	\$49.90
	1031B	500 ml	\$92.70
Chromotrope 2R Stain, Aqueous (Use: Mod. of Wheatley Trichrome Stain for microsporidia, etc.)	10351A	250 ml	\$80.20
Congo Red Stain Set, Bennhold, Amyloid Solution A: Congo Red Stain 1%, Aqueous Solution B: Alkaline Alcohol	1038A	500 ml set	\$80.40 H
Congo Red Stain Set, Puchtler, Amyloid (Use: Mod. of Bennhold; reduces background birefringence) Solution A: Sodium Hydroxide 1%, Aqueous, 25 ml or 50 ml Solution B: Congo Red Stain, Alcoholic, 250 ml or 500 ml	1037A	250 ml set	\$70.80 H
	1037B	500 ml set	\$127.80 H
Crystal Violet, Powder	1040A	10 g	\$13.50
	1040B	25 g	\$23.00
	1040C	100 g	\$56.00
Crystal Violet Stain 1%, Aqueous, Brown-Hopps (Use: Brown-Hopps mod. Gram Stain)	1041A	500 ml	\$34.30
Crystal Violet-Oxalate Stain, Alcoholic, Brown-Brenn (Use: Brown-Brenn Gram Stain)	1042A	500 ml	\$26.40
Crystal Violet-Oxalate Stain, Alcoholic, Hucker-Twort (Use: Hucker-Twort Gram Stain)	10422A	250 ml	\$34.70
	10422B	500 ml	\$47.80
Crystal Violet Stain, Lieb, Alcoholic Easy, backup Amyloid Stain	10421A	100 ml	\$61.70
Eosin Stain 0.5%, Aqueous NaCl (Use: Sperm Vitality Stain)	1068B	500 ml	\$39.60



PRODUCT NAME	PART	VOLUME	PRICE	
Fast Green Stain 2.5%, Aqueous	10852A	250 ml	\$50.40	
	10852B	500 ml	\$73.40	
Fouchet Reagent (Use: Hall Method for Bilirubin)	1095A	250 ml	\$76.70	
	1095B	500 ml	\$124.40	
Giemsa, Powder	1110A	10 g	\$36.30	
	1110B	25 g	\$87.50	
Giemsa Stock Stain Dilute with water or buffer.	1120A	500 ml	\$62.50	H
	1120B	1 L	\$79.80	H
Giemsa Stock Stain, Romanowski (Use: The traditional Giemsa for Hematology)	11215A	500 ml	\$37.70	H
	11215A	6 x 500 ml	\$129.60	H
Giemsa Stock Stain, Wolbach (Use: Wolbach & May-Grunwald mod. Giemsa Stain)	1121A	500 ml	\$52.10	H
	1121B	1 L	\$90.30	H
Gomori Basement Membrane Silver Set Solution A: Methenamine 3%, Aqueous, 250 ml Solution B: Silver Nitrate 5%, Aqueous, 50 ml Solution C: Sodium Borate 5%, Aqueous, 50 ml	1131A	250 ml set	\$73.70	
Grocott Methenamine Silver Set, GMS (Use: Fungus Stain) Shelf Life: 9 mos. Store at 4°C. Solution A: Silver Nitrate Solution B: Methenamine Borate Each bottle 500 ml Each bottle 1 L	1142A	1L set	\$97.20	
	1142B	2L set	\$176.20	
	1162A	250 ml	\$35.00	H
	1162B	500 ml	\$65.50	
Hematoxylin 3%, Alcoholic (Use: Elastic-Masson Trichrome combination stain with myocardial applications)	11623B	500 ml	\$158.50	H
	11623C	1 L	\$221.20	H
Hematoxylin 5%, Alcoholic (Use: Verhoeff-Van Gieson Elastic Stain)	1210A	500 ml	\$32.50	H
	1210B	1 L	\$58.40	H
Jenner Stock Stain (Use: May-Grunwald mod. Giemsa Stain)	12204A	250 ml	\$21.20	
	12204B	500 ml	\$28.00	
Light Green SF Yellowish Stain 0.02%, Aqueous (Use: Counterstain for GMS Fungus Stain and McLetchie Light Green Trichrome)	12203A	250 ml	\$20.70	
	12203B	500 ml	\$28.40	
Light Green SF Yellowish Stain 0.2%, Aqueous (Use: Stock solution for counterstain for GMS fungus)	12202A	250 ml	\$23.00	
	12202B	500 ml	\$30.40	
Light Green SF Yellowish Stain 1%, Aqueous (Use: Stock solution for GMS & Gomori Basement Membrane)	1220C	250 ml	\$31.80	
	1220B	500 ml	\$59.10	

STAINS & STAIN KITS

PRODUCT NAME	PART	VOLUME	PRICE	
Special Stains				
Light Green SF Yellowish Stain 2%, Aqueous (Use: Elastic-Masson Trichrome and stock solutions for GMS & Gomori Basement Membrane)	1221A	250 ml	\$48.70	
	1221B	500 ml	\$71.20	
Luxol Fast Blue Stain Set (Use: Myelin Stain) Solution A: Luxol Fast Blue Stain 0.1%, Alcoholic Solution B: Lithium Carbonate 1%, Aqueous	12218A	500 ml set	\$70.70	H
	12218B	1 L set	\$100.80	H
Metanil Yellow Stain, Aqueous (Use: Counterstain for Mucicarmine)	12235B	250 ml	\$39.80	
	12235C	500 ml	\$72.20	
Methyl Green Stain 2%, Aqueous (Use: Counterstain for IHC procedures; Not chloroform extracted)	12253A	250 ml	\$81.70	
	12253B	500 ml	\$150.40	
	12253C	1 L	\$232.00	
Methylene Blue, Powder	1230A	10 g	\$13.20	
	1230B	25 g	\$15.70	
	1230C	100 g	\$39.10	
Methylene Blue Stain 0.5%, Working Aqueous (Use: AFB Counterstain for Ziehl-Neelson, Fite)	12402A	250 ml	\$32.00	
Methylene Blue Stain 1.4%, Alcoholic (Use: AFB Counterstain)	1240C	250 ml	\$27.20	H
	1240A	500 ml	\$30.10	H
	1240B	1 L	\$36.60	H
Mucicarmine Stock Stain, Mayer Shelf Life: 6 mos. Store at 4°C.	1250A	125 ml	\$59.70	H
	1250B	500 ml	\$114.60	H
	1250C	1 L	\$155.50	H
New Methylene Blue N Stain, Aqueous (Use: Reticulocyte Count for Hematology)	1270B	500 ml	\$29.50	
	1270C	1 L	\$40.40	
Nigrosin, Powder	1271A	25 g	\$12.10	
Nigrosin Stain 10%, Aqueous (Use: Mod. of sperm vitality stain)	1272B	500 ml	\$90.30	
Nuclear Fast Red Stain, Kernechtrot (Use: Counterstain for Iron Stains, Schmorl Melanin & IHC Stains)	1255A	250 ml	\$61.70	
	1255C	500 ml	\$91.30	
	1255B	1 L	\$120.10	
Oil Red O, Powder	1276B	25 g	\$19.30	
Oil Red O Stain, Propylene Glycol (Use: Stains neutral lipids)	12772A	250 ml	\$43.50	
	12772B	500 ml	\$73.90	
Oil Red O Stain, Isopropanol (Use: Stains neutral lipids)	1277A	250 ml	\$48.70	H
	1277B	500 ml	\$65.10	H



PRODUCT NAME	PART	VOLUME	PRICE	
Papanicolaou Stain, EA-36, Fast Green Modified formula with Fast Green.	1300A	500 ml	\$20.20	H
	1300B	1 gal	\$50.20	H
	1300C	10 L	\$94.80	H
Papanicolaou Stain, EA-36, Light Green Original formula with Light Green.	1302A	500 ml	\$19.00	H
	1302B	1 gal	\$52.30	H
	1302C	10 L	\$125.00	H
Papanicolaou Stain, EA-50, Fast Green Modified formula with Fast Green.	1310A	500 ml	\$21.10	H
	1310B	1 gal	\$52.30	H
	1310C	10 L	\$97.30	H
Papanicolaou Stain, EA-50, Light Green Original formula with Light Green.	1312A	500 ml	\$21.10	H
	1312B	1 gal	\$52.30	H
	1312C	10 L	\$93.00	H
Papanicolaou Stain, EA-65, Fast Green Modified formula with Fast Green.	1320A	500 ml	\$20.20	H
	1320B	1 gal	\$50.20	H
	1320C	10 L	\$97.20	H
Papanicolaou Stain, EA-65, Light Green Original formula with Light Green.	1322A	500 ml	\$19.00	H
	1322B	1 gal	\$47.10	H
	1322C	10 L	\$101.30	H
Papanicolaou Stain, OG-6	1330A	500 ml	\$21.10	H
	1330B	1 gal	\$53.20	H
	1330C	10 L	\$101.30	H
Phosphotungstic Acid Hematoxylin (PTAH) Stain (Use: Demonstrates muscle striations)	1334A	500 ml	\$67.50	
	1334B	1 L	\$120.10	
Rhodamine B, Powder (Use: IF Stain for AFB)	13397A	10 g	\$18.20	
Rhodanine, Powder (Use: Copper Stain)	1053A	5 g	\$43.00	
Rhodanine Stock Stain 0.2%, Alcoholic (Use: Stock solution for Copper Stain) Store at 4°C.	10531A	250 ml	\$41.50	H
	10531B	500 ml	\$74.60	H
Saffron, Powder	1362A	5 g	\$293.80	
Safranin O, Powder	1340A	10 g	\$13.30	
	1340B	25 g	\$23.50	
	1340C	100 g	\$51.00	
Safranin O Stain 0.25%, Aqueous (Use: Counterstain for some Gram methods)	1360A	500 ml	\$23.20	
	1360C	1 L	\$30.00	
	1360E	1 gal	\$72.00	
Safranin O Stain 1%, Aqueous (Use: Cartilage Staining)	1350A	125 ml	\$10.10	
	1350B	500 ml	\$36.30	

STAINS & STAIN KITS

PRODUCT NAME	PART	VOLUME	PRICE
Special Stains			
Schiff Reagent, McManus Store at room temperature for up to one year.	1371A	125 ml	\$32.40
	1371B	500 ml	\$104.90
	1371C	1 L	\$124.70
	1371D	1 gal	\$256.30
Silver Nitrate, ACS Powder	13803A	25 g	\$129.60
	13803B	100 g	\$432.00
Silver Nitrate 1%, Aqueous (Use: Steiner-Steiner & Steiner-Chapman Modified Silver Stain)	13804A	250 ml	\$66.30
	13804B	500 ml	\$124.70
Silver Nitrate 5%, Aqueous (Use: Von Kossa Calcium Stain, Gomori Basement Membrane and Gomori Methenamine Silver for Urates)	13805A	250 ml	\$119.30
	13805B	500 ml	\$157.80
Silver Nitrate 10%, Aqueous (Use: Fontana-Masson Silver Stain for argentaffin granules, Gordon & Sweets Reticulum)	13806A	250 ml	\$165.80
	13806B	500 ml	\$273.80
Silver Nitrate 20%, Aqueous	13807	250 ml	\$273.00
Sudan Black B Stain, Propylene Glycol (Use: Stains Neutral Lipids & Phospholipids)	1401A	125 ml	\$22.60
	1401B	250 ml	\$42.10
Sudan III Stain, Saturated Alcoholic	1390A	500 ml	\$37.70 H
Sudan IV Stain, Herxheimer Alcoholic	1400A	125 ml	\$30.60 H
Tartrazine Stain 0.25%, Acetic Aqueous (Use: Counterstain for Mucicarmine)	14016A	250 ml	\$37.60
	14016B	500 ml	\$59.40
Tartrazine Stain 1.5%, Aqueous (Use: Brown-Hopps Gram Stain)	14015A	500 ml	\$21.80
	14015B	1 L	\$31.30
Thiazine Stain	10522A	500 ml	\$26.50
	10522B	1 gal	\$83.20
Toluidine Blue Stain 0.1%, Aqueous (Use: Mast cell stain and MOHS basal cell carcinoma)	14027A	250 ml	\$40.70
	14027B	500 ml	\$68.80
	14027D	1 gal	\$295.10
Toluidine Blue Stock Stain 1%, Aqueous (Use: Toluidine Blue/Alcian Yellow Helicobacter Stain)	14028A	500 ml	\$47.70
	14028B	1 L	\$83.70
Trichrome Stain, Gomori One-Step, Aniline Blue Store at 4°C.	1403C	250 ml	\$42.60
	1403B	500 ml	\$73.50
Trichrome Stain, Gomori One-Step, Fast Green Store at 4°C.	14021A	250 ml	\$49.90
	14021B	500 ml	\$78.00
Trichrome Stain, Gomori One-Step, Light Green Store at 4°C.	1402C	250 ml	\$46.20
	1402B	500 ml	\$73.40



PRODUCT NAME	PART	VOLUME	PRICE	
Twort Stain Set (Use: Hucker-Twort Gram Stain) Solution A: Neutral Red Stain 1%, Alcoholic, 250 ml or 500 ml Solution B: Fast Green Stain 1%, Alcoholic, 100 ml or 200 ml	14034A 14034B	250 ml set 500 ml set	\$56.00 \$93.20	H H
Van Gieson Stain (Use: Counterstain for VVG Elastic Stain)	1404A 1404B	500 ml 1 L	\$63.70 \$106.40	
Victoria Blue Stain, Alcoholic (Uses include: Hepatitis B surface antigen, copper and venous invasion)	1406A 1406C	500 ml 1 L	\$142.40 \$258.80	H H
Wright, Powder	1410A 1410B 1410C	10 g 25 g 100 g	\$17.20 \$25.50 \$84.60	
Wright Stain (Use: Peripheral and BM smears)	1420A 1420B 1420C	500 ml 1 L 1 gal	\$26.40 \$29.60 \$78.80	H H H
Wright Stain, Buffered (Use: Peripheral and BM smears) Faster staining procedure with more polychromasia than unbuffered.	1422A 1422B 1422C	500 ml 1 L 1 gal	\$25.40 \$30.90 \$78.80	H H H
Wright Stain, Modified (Use: Stronger for auto stainers)	1421A 1421B 1421C	500 ml 1 L 1 gal	\$24.80 \$40.60 \$78.80	H H H
Xanthene Stain	10521A 10521B	500 ml 1 gal	\$34.70 \$72.30	

Stain Kits

Fat, Oil Red O, Propylene Glycol Stain Kit

9119A 250 ml kit \$186.50

Carbohydrate Stains

Alcian Blue 1%, pH 2.5 Stain Kit Shelf Life: 6 mos.

9102A 250 ml kit \$149.70
9102B 500 ml kit \$191.90

Alcian Blue/PAS Stain Kit Shelf Life: 6 mos.

91022A 250 ml kit \$179.10
91022B 500 ml kit \$323.20

Colloidal Iron, Müller-Mowry Stain Kit

9110A 250 ml kit \$272.70 H

Mucin, Mayer Mucicarmine Stain Kit Shelf Life: 6 mos.

9151A 250 ml kit \$147.00 H
9151B 500 ml kit \$222.50 H

Periodic Acid Schiff (PAS) Stain Kit

9162A 250 ml kit \$196.10 H
9162B 500 ml kit \$272.40 H

STAINS & STAIN KITS

PRODUCT NAME	PART	VOLUME	PRICE	
Connective Tissue & Muscle Fiber Stains				
Basement Membrane, Gomori Stain Kit	9167A	250 ml kit	\$312.80	
	9167B	500 ml kit	\$386.30	
Phosphotungstic Acid Hematoxylin (PTAH) Stain Kit	9111A	250 ml kit	\$137.20	H
Elastic, Verhoeff Stain Kit	9116A	250 ml kit	\$190.60	H
	9116B	500 ml kit	\$257.50	H
Movat-Russell Modified Pentachrome Stain Kit	9150A	250 ml kit	\$467.00	H
	9150B	500 ml kit	\$752.50	H
Reticulum, Gordon & Sweets Stain Kit	9168A	250 ml kit	\$376.30	
	9168B	500 ml kit	\$710.20	H
Trichrome, Gomori One-Step, Aniline Blue Stain Kit	9176B	250 ml kit	\$138.90	H
	9176A	500 ml kit	\$212.80	H
Trichrome, Gomori One-Step, Fast Green Stain Kit	9175A	250 ml kit	\$138.90	H
	9175B	500 ml kit	\$212.80	H
Trichrome, Masson, Aniline Blue Stain Kit	9179A	250 ml kit	\$294.70	H
	9179B	500 ml kit	\$381.20	H
Trichrome, Masson, Fast Green Stain Kit	9180A	250 ml kit	\$273.10	H
	9180B	500 ml kit	\$462.70	H
Trichrome, McLetchie, Aniline Blue Stain Kit	9177A	250 ml kit	\$182.00	H
	9177B	500 ml kit	\$302.10	H
Trichrome, McLetchie, Light Green Stain Kit	9178A	250 ml kit	\$182.00	H
	9178B	500 ml kit	\$302.10	H
Microorganism Stains				
AFB, Basic Fuchsin, Ellis & Zabrowarny Stain Kit	91012A	250 ml kit	\$102.80	H
	91012B	500 ml kit	\$136.90	H
AFB, Fite Stain Kit	91013A	250 ml kit	\$134.80	H
AFB, Ziehl-Neelsen Stain Kit	9101A	250 ml kit	\$67.90	H
Differential Stain Kit (Use: Hematological Staining & Special Stain for <i>Helicobacter pylori</i>)	9112B	500 ml kit	\$68.50	H
Fungus, Grocott Methenamine Silver (GMS) Stain Kit Shelf Life: 9 mos.	9121A	250 ml kit	\$249.90	H
	9121B	500 ml kit	\$455.30	H
Fungus, PAS/Light Green Stain Kit	9122A	250 ml kit	\$145.40	
	9122B	500 ml kit	\$267.30	
Gram, Brown-Hopps Stain Kit	9124A	250 ml kit	\$132.00	H
Gram, Brown-Brenn, Modified Stain Kit	9123A	250 ml kit	\$147.80	H



PRODUCT NAME	PART	VOLUME	PRICE	
Gram, Hucker-Twort Stain Kit	9125A	250 ml kit	\$147.90	H
Helicobacter, Toluidine Blue / Alcian Yellow Stain Kit	9130A	250 ml kit	\$244.60	H
	9130B	500 ml kit	\$363.60	H
Steiner-Chapman Modified Silver Stain Kit No Uranyl Nitrate!	9172A	250 ml kit	\$219.20	H
	9172B	500 ml kit	\$386.30	H
Steiner-Steiner Modified Silver Stain Kit	9171A	250 ml kit	\$430.50	H
	9171B	500 ml kit	\$568.60	H
Neurological Stains				
Bielschowsky, Lester-King Modified Stain Kit	9154A	250 ml kit	\$314.50	H
	9154B	500 ml kit	\$572.30	H
Luxol Fast Blue (LFB) - Cresyl Violet Stain Kit	9155A	250 ml kit	\$224.10	H
	9155B	500 ml kit	\$395.50	H
Pigment, Substance & Mineral Stains				
Amyloid, Bennhold Congo Red Stain Kit	9103A	250 ml kit	\$95.40	H
Amyloid, Puchtler Congo Red Stain Kit	9104A	250 ml kit	\$128.80	H
Copper, Rhodanine Stain Kit	9113A	250 ml kit	\$109.40	H
Fontana Masson Stain Kit	9105A	250 ml kit	\$422.80	H
Iron, Gomori Prussian Blue Stain Kit	9136A	250 ml kit	\$121.60	H
	9136B	500 ml kit	\$188.00	H
Buffers				
Hyaluronidase Buffer (Use: Buffer for digestion in Alcian Blue pH 2.5; Does not include enzyme) Shelf Life: 6 mos. Store at 4°C.	1150B	1 L	\$52.80	
	1150C	1 gal	\$90.30	
Phosphate Buffer, pH 6.0 (Use: Buffer for Glycogen Digestion; Does not include enzyme) Shelf Life: 6 mos. Store at 4°C.	13312A	500 ml	\$31.60	
	13312B	1 L	\$52.50	
	13312C	1 gal	\$111.50	
Phosphate Buffer, pH 7.0 (Use: With Giemsa Stain) Shelf Life: 6 mos. Store at 4°C.	1331A	500 ml	\$16.90	
	1331B	1 L	\$22.50	
Wright Stain Buffer, pH 6.8 Shelf Life: 6 mos. Store at 4°C.	1430A	500 ml	\$24.70	
	1430B	1 L	\$37.20	
	1430C	1 gal	\$68.30	

CHEMICALS & REAGENTS

PRODUCT INFORMATION

- Shelf Life is **one year** from date manufactured unless otherwise stated.
- Store at 15-30°C unless otherwise stated.
- Acetone, Alcohols and Xylene do not have expiration dates. Shelf Life is 2 yrs. once opened.

H = Hazardous

PRODUCT NAME

PART

VOLUME

PRICE

Alcohols

Alcohol, Reagent ACS, 100%	10051A	500 ml	\$14.70	H
Ethyl Alcohol 90%	10051A	6 x 500 ml	\$49.80	H
Methanol Alcohol 5%	10051C	1 gal	\$35.40	H
Isopropyl Alcohol 5%	10051C	4 x 1 gal	\$70.80	H
Alcohol, Reagent ACS, 95%	10052C	1 gal	\$35.40	H
	10052C	4 x 1 gal	\$70.80	H
Alcohol, Reagent ACS, 80%	10053C	1 gal	\$35.40	H
	10053C	4 x 1 gal	\$70.80	H
Alcohol, Reagent ACS, 70%	10054C	1 gal	\$35.40	H
	10054C	4 x 1 gal	\$70.80	H
Alcohol, Ethyl Denatured, 100%	10841A	1 gal	\$41.30	H
Ethyl Alcohol 86%	10841B	1 gal	\$70.80	H
Methyl Alcohol 4%	10841C	5 gal can	\$81.60	H
Isopropyl Alcohol 9%				
Methyl Isobutyl Ketone 1%				
Alcohol, Ethyl Denatured, 95%	10842A	1 gal	\$34.70	H
	10842A	4 x 1 gal	\$70.80	H
	10842C	5 gal can	\$81.60	H
Alcohol, Ethyl Denatured, 80%	10843A	1 gal	\$34.70	H
	10843A	4 x 1 gal	\$70.80	H
	10843C	20 L cube	\$81.60	H
Alcohol, Ethyl Denatured, 70%	10844A	1 gal	\$34.70	H
	10844A	4 x 1 gal	\$70.80	H
	10844C	20 L cube	\$81.60	H
Alcohol, Isopropyl ACS, 100%	12094A	500 ml	\$15.80	H
	12094A	6 x 500 ml	\$58.80	H
	12094E	4 L	\$58.80	H
	12094E	4 X 4 L/cs	\$188.40	H
	12094F	5 gal can	Inquire	H
Alcohol Isopropyl, 70%	12093C	4 x 1 gal	\$182.80	H



PRODUCT NAME	PART	VOLUME	PRICE	
Alcohol, Methanol Anhydrous, ACS (Use: Fixative for blood smears)	12236A	500 ml	\$15.70	H
	12236A	6 x 500 ml	\$45.60	
	12236C	1 gal	\$45.30	H
	12236C	4 x 1 gal	\$122.80	
Eosin Y 1%, Isopropyl Alcohol (Use: In Leica Peloris Processor as tint for tissue prior to embedding)	1071A	250 ml	\$37.40	H
	1071B	500 ml	\$41.70	H
	1071C	1 L	\$57.20	H

Decalcifying Solutions

Decalcifying Solution, EDTA/Sucrose	1048B	1 L	\$51.70	
	1048C	1 gal	\$182.30	
	1048D	10 L cube	\$326.90	
Decalcifying Solution, Formic Acid 5%, Aqueous	1049B	1 L	\$28.10	
	1049C	1 gal	\$100.10	
	1049E	20 L cube	\$210.70	
Decalcifying Solution, HCL/EDTA (Use: Excellent for bone marrows) Shelf Life: 2 yrs.	1050A	1 L	\$21.40	
	1050A	6 x 1 L	\$101.40	
	1050C	1 gal	\$63.90	
Decalcifying Solution, Formic/Citrate Shelf Life: 2 yrs.	10492A	500 ml	\$29.50	H
	10492B	1 L	\$35.60	H
	10492C	1 gal	\$89.00	H
Decalcifying Solution, Formic Acid/Formalin Shelf Life: 2 yrs.	10493B	1 L	\$24.20	
	10493C	1 gal	\$66.30	
Decalcification End Point Set (Use: Determination of complete decalcification of specimen) Solution A: Ammonium Hydroxide 5%, Aqueous Solution B: Ammonium Oxalate 5%, Aqueous	1051A	250 ml each bottle	\$43.90	

CHEMICALS & REAGENTS

PRODUCT NAME	PART	VOLUME	PRICE
Fixatives			
B-5 Mercury Substitutes			
AZF Fixative (Use: Replacement for mercury based fixative) Working Solution. Mercury Free!	1009A	1 L	\$23.20
	1009B	1 gal	\$76.40
	1009C	20 L cube	\$104.10
AZF Fixative Vial 40 ml vial, 20 ml fill 60 ml vial, 30 ml fill	10091D	100/cs	\$132.80
	10091C	100/cs	\$137.90
B-5 Fixative Modified, Zinc Chloride (Use: Alternative to mercury based fixative) Mercury free!	1015A	1 L	\$47.30
	1015A	6 X 1 L	\$249.00
	1015C	1 gal	\$93.70
Electron Microscopy			
Carson Modified Millonig Formalin (Use: Light and EM Fixative) 30 ml vial, 15 ml fill	12445A	1 gal	\$10.10
	12445C	20 L cube	\$49.10
	12445E	100/cs	\$62.50
Glutaraldehyde 50%, Reagent Grade (Use: Electron Microscopy)	1123A	500 ml	\$108.60 H
Millonig Buffer, pH 7.0 (Use: Electron Microscopy)	12442A	500 ml	\$72.20
	12442B	1 L	\$105.70
Zamboni Fixative (Use: For Light and Electron Microscopy)	1459A	1 L	\$92.20
Formalin			
Formalin 10%, Phosphate Buffered	1090N	1 gal	\$11.00
	1090N	4 x 1 gal	\$35.20
	1090K	20 L cube	\$37.10
Formalin Concentrate Form-C (Buffered) One gallon of concentrate makes 5 gallons of 10% NBF. Thoroughly mixes in 30 seconds.	ABFORM-04	4 x 1 gal/cs	\$118.00
Formaldehyde 37-40%, ACS	1089A	500 ml	\$19.30 H
	1089C	1 gal	\$45.20 H
	1089B	20 L cube	\$84.90 H
Formalin, Alcoholic, 10/90	10932A	4 x 1 gal	\$86.80 H
	10932B	20 L cube	\$86.80 H
Formalin, Alcoholic, 20/80	10933A	4 x 1 gal	\$75.80 H
	10933B	20 L cube	\$90.30 H



PRODUCT NAME	PART	VOLUME	PRICE	
Formalin 23%, Phosphate Buffered (Use: CNS tissue e.g., large animal brains)	1091A	1 gal	\$24.30	H
	1091C	10 L cube	\$59.10	H
Lymph Node Grossing Aid (Formaldehyde in Ethanol, Acetic Acid) (Use: Enhances gross visibility of lymph nodes in large resections)	1093A	4 x 1 gal	\$126.20	H
	1093B	20 L cube	\$133.20	H
Formalin Vial 10%, Phosphate Buffered				
20 ml vial, 10 ml fill	1090A	100/cs	\$27.60	
40 ml vial, 20 ml fill	1090B	100/cs	\$36.40	
32 oz vial, 500 ml fill	1090E	9/cs	\$23.10	
2 L vial, 1 L fill	1090F	4/cs	\$16.10	
FrozFix®				
FrozFix® Starter Set (Each bottle 250 ml) Solution A: Reagent Alcohol 95% Solution B: Acetone-Alcohol Solution C: FrozFix® Fixative Solution D: Tris Buffered Saline + Tween, Mod. 0.05M, pH7.6, 10X	10972	250 ml set	\$44.80	H
FrozFix® Fixative	1096A	250 ml	\$16.30	H
	1096B	500 ml	\$31.90	H
	1096C	1 L	\$44.70	H
	1096D	1 gal	\$118.90	H
	1096E	20 L cube	Inquire	H
FrozFix® Conditioning Set (Necessary step prior to FrozFix® Fixative) Solution A: Reagent Alcohol 95% Solution B: Acetone-Alcohol	1097A	250 ml set	\$11.30	H
	1097B	500 ml set	\$13.70	H
	1097C	1 L set	\$18.50	H
	1097D	1 gal set	\$59.50	H
Tris Buffered Saline + Tween Mod. 0.05M, pH 7.6, 10X (Use: FrozFix® IHC Applications) Shelf Life: 6 mos. Store at 4° C.	140305A	500 ml	\$19.70	
	140305B	1 L	\$32.50	
	140305C	1 gal	\$98.20	
	140305D	10 L cube	Inquire	
	140305E	20 L cube	Inquire	
Tris Buffered Saline Mod. 0.05M, pH 7.6, 10X (Use: FrozFix® IHC Applications)	140306A	500 ml	\$18.50	
	140306B	1 L	\$34.90	
	140306C	1 gal	\$96.10	
	140306D	10 L cube	Inquire	
	140306E	20 L cube	Inquire	

CHEMICALS & REAGENTS

PRODUCT NAME	PART	VOLUME	PRICE
Michel Transport Medium			
Michel Transport Medium (Use: Preserves morphology and antigenicity during transport of tissues) Shelf Life: 8 mos. Store at 4° C.	1242A	500 ml	\$47.30
	1242B	1 L	\$65.10
	1242C	1 gal	\$149.80
Michel Transport Medium Vial Shelf Life: 8 mos. Store at 4° C.			
30 ml vial, 15 ml fill	12423C	50/cs	\$77.60
30 ml vial, 15 ml fill	12423A	100/cs	\$144.70
60 ml vial, 30 ml fill	12423B	100/cs	\$152.00
Michel Wash Solution Shelf Life: 8 mos. Store at 4° C.	1243A	500 ml	\$38.80
	1243B	1 L	\$49.60
	1243C	1 gal	\$131.90
Other Fixatives			
Bouin Fluid	1020A	1 L	\$36.50
	1020B	1 gal	\$74.10
	1020C	20 L cube	\$149.50
Bouin Fluid Vial Ship Ground Only.			
40 ml vial, 20 ml fill	10201B	100/cs	\$144.70
60 ml vial, 20 ml fill	10201E	72/cs	\$107.60
60 ml vial, 40 ml fill	10201C	72/cs	\$114.60
Davidson Fixative	1045A	1 gal	\$100.30 H
Differential Stain Fixative Ship Ground Only.	10523A	500 ml	\$31.40 H
	10523B	1 gal	\$76.70 H
Formic Acid 96%, ACS	10943A	250 ml	\$33.40 H
	10943B	500 ml	\$64.20 H
	10943C	2 L	\$171.50 H
Hollande Fixative (Use: Recommended for GI Biopsies)	1208A	1 L	\$57.90
	1208B	1 gal	\$133.30
Hollande Fixative Vial Ship Ground Only.			
40 ml vial, 20 ml fill	12082A	100/cs	\$144.70
Karnovsky Fixative (Use: EM; Rapid penetration, esp. neurological tissues) Shelf Life: 6 mos. Store at 4° C.	1215A	1 L	\$117.90
	1215B	1 gal	\$223.00
Karnovsky Fixative, Phosphate Buffered Shelf Life: 6 mos. Store at 4° C.	12158A	1 L	\$133.50



PRODUCT NAME	PART	VOLUME	PRICE	
Karnovsky Fixative Vial Shelf Life: 6 mos. Store at 4° C.				
30 ml vial, 15 ml fill	12151A	50/cs	\$263.30	
Newcomer Fixative (Use: Mucopolysaccharides and nucleoproteins)	1253A	1 L	\$113.00	
Paraformaldehyde 10%, Buffered (Use: Cell flow fixation)	13301A	500 ml	\$87.80	
	13301B	1 L	\$108.90	
	13301C	1 gal	\$262.70	
Phenol-Formalin Fixative	13307A	1 gal	\$106.40	H
	13307A	4 X 1 gal	\$389.60	H
Saccomano Collection Fluid (Use: For collection and fixation of mucoid specimens)	1361B	1 gal	\$30.90	H
60 ml vial, 30 ml fill	1361D	50/cs	\$86.80	H
Zamboni Fixative, Modified (Use: Light and EM fixative)	13381A	1 L	\$66.00	
	13381B	1 gal	\$123.70	
Zenker Fixative, Modified, Zinc Chloride	1461A	1 L	\$43.00	
	1461B	1 gal	\$107.00	
Zinc Formalin Fixative, Unbuffered Zinc Sulfate (Use: On auto processor or as a prefixative) Working solution.	1482A	1 L	\$13.20	
	1482B	1 gal	\$22.40	
	1482C	20 L cube	\$86.80	

Xylenes and Substitutes

Xylene, Histology Grade	1446A	500 ml	\$16.70	H
	1446C	1 gal	\$36.90	H
	1446C	4 x 1 gal	\$88.80	H
	1446E	5 gal can	\$97.90	H
	1446F	55 gal	Inquire	H
Xylene, ACS	1445A	500 ml	\$19.30	H
	1445A	6 x 500 ml	\$60.00	H
	1445C	1 gal	\$61.10	H
	1445C	4 x 1 gal	\$146.80	H
	1445E	5 gal can	\$104.10	H
	1445F	55 gal	Inquire	H
Xylene / Peanut Oil, 2:1 (Use: Fite Stain for Leprosy and Nocardia)	1449A	250 ml	\$39.80	H
	1449B	500 ml	\$55.50	H
	1449C	1 L	\$91.80	H
Slide Brite™ An extremely safe, odorless alternative to Xylene.	ABSB-04	4 x 1 gal/cs	\$164.00	

CHEMICALS & REAGENTS

PRODUCT NAME	PART	VOLUME	PRICE	
Chemicals				
Acetic Acid 0.5%, Aqueous (Use: Trichrome, Gomori One-Step)	100121A	500 ml	\$10.00	
	100121B	1 L	\$21.30	
Acetic Acid 1%, Aqueous (Use: Masson and Gomori Trichrome and May-Grunwald Giemsa)	10012A	500 ml	\$17.50	
Acetic Acid 3%, Aqueous (Use: Alcian Blue, pH 2.5)	10017A	250 ml	\$13.80	
	10017B	500 ml	\$16.00	
Acetic Acid, Glacial, ACS	10010A	500 ml	\$29.90	H
	10010A	6 x 500 ml	\$90.90	H
Acetone, ACS (Use: Various, including Brown-Brenn and Hucker-Twort Gram Stains)	10014A	500 ml	\$15.50	H
	10014A	6 x 500 ml	\$58.20	
	10014B	1 L	\$29.60	H
	10014D	4 L	\$53.80	H
	10014D	4 x 4 L	\$148.20	H
	10014F	5 gal can	\$115.20	H
Acetone-Alcohol 1:1 (Use: Gram Stain Brown-Brenn)	10016A	250 ml	\$16.40	H
	10016B	500 ml	\$22.60	H
Acetone-Xylene 1:1 (Use: Brown-Brenn and Brown-Hopps Gram Stains)	10015A	500 ml	\$15.30	H
Acid-Alcohol 1% (Use: H&E, AFB, PAS and Fite Stains)	10011A	500 ml	\$16.10	H
	10011B	1 L	\$17.10	H
	10011C	1 gal	\$45.40	H
	10011C	4 x 1 gal	\$157.20	
Acidulated Water pH, 4.0-4.1 (Use: Warthin-Starry Method)	10013A	500 ml	\$20.00	
	10013B	1 L	\$33.10	
Ammonium Hydroxide 28-30%, ACS (Use: Reticulum, Gordon and Sweets; Lester King, Bielschowsky)	1006A	500 ml	\$38.00	H
	1006A	6 x 500 ml	\$78.60	H
Borax (Sodium Borate) 5%, Saturated Alcoholic (Use: Sulfated Alcian Blue for Amyloid)	1019A	500 ml	\$34.10	H
	1019B	1 L	\$59.50	H
Chromic Acid 10%, Aqueous (Use: Stock solution for Grocott Solution GMS)	1034A	500 ml	\$75.90	H
Chromic Acid 5%, Aqueous (Use: Working dilution for Grocott Solution GMS)	10341A	250 ml	\$32.50	H
	10341B	500 ml	\$63.70	H
Colloidal Iron Stock (Use: Mowry Colloidal Iron)	10365B	500 ml	\$108.80	
	10365C	1 L	\$217.60	



PRODUCT NAME	PART	VOLUME	PRICE	
Ferric Chloride 1%, Aqueous (Use: Schmorl Melanin Stain)	10855A	250 ml	\$16.10	
	10855B	500 ml	\$23.70	
Ferric Chloride 10%, Aqueous (Use: Verhoeff Elastic and Hall Bile Stains)	10856A	250 ml	\$28.10	
	10856B	500 ml	\$39.40	
Ferric Chloride 2%, Aqueous (Use: Working solution for VVG Elastic Stain)	108553A	250 ml	\$16.90	
	108553B	500 ml	\$28.60	
Ferric Chloride 29%, Aqueous Stock (Use: Weigert, Lillie and Mallory Iron Hematoxylin)	1086A	500 ml	\$38.30	
	1086B	1 L	\$70.70	
Gallego Solution (Use: Brown-Hopps modified Gram Stain)	1098A	500 ml	\$24.00	
	1098B	1 L	\$42.50	
Gold Chloride 0.1%, Aqueous (Use: GMS Stain)	11285A	100 ml	\$99.40	
	11285B	250 ml	\$131.00	
Gold Chloride 0.2%, Aqueous (Use: Gordon and Sweets Reticulum, Fontana-Masson Argentaffin and GMS Stains)	11286A	100 ml	\$109.30	
	11286B	250 ml	\$149.90	
Gold Chloride 0.25%, Aqueous (Use: Gomori Basement Membrane Stain)	11287A	100 ml	\$115.80	
	11287B	250 ml	\$153.50	
Gold Chloride 1%, Aqueous (Use: Gordon and Sweets and Gomori Basement Membrane, Fontana-Masson Argentaffin and GMS Stains)	1129A	100 ml	\$171.10	
Gum Mastic 2.5%, Alcoholic	1145A	250 ml	\$39.40	H
	1145B	500 ml	\$60.50	H
Hydrochloric Acid 0.1 N, Aqueous (Use: Alcian Blue 1.0 pH Sulfated Mucosubstances)	1207A	250 ml	\$10.10	
	1207B	500 ml	\$14.00	
Hydrochloric Acid 20%, Aqueous (Use: Gomori Mod. Iron Stain, Prussian Blue and stock for Colloidal Iron)	12087A	500 ml	\$28.90	H
	12087B	1 L	\$36.50	H
Hydrochloric Acid 36-38%, ACS	12084A	500 ml	\$34.30	H
	12084A	6 x 500 ml	\$165.00	H
	12084C	1 gal	\$51.40	H
Hydrochloric Acid 5%, Aqueous (Use: For acid cleaning glassware prior to silver, calcium and iron staining)	12086A	1 L	\$24.70	H
	12086B	1 gal	\$63.70	H
	12086C	20 L	\$181.00	H
Hydrogen Peroxide 30%, Aqueous (Use: IHC procedures)	1206A	16 oz	\$53.30	H

CHEMICALS & REAGENTS

PRODUCT NAME	PART	VOLUME	PRICE	
Chemicals				
Iodine, Gram, Aqueous (Use: Brown-Brenn and Brown-Hopps Gram Stains and mercury pigment removal)	1140A	500 ml	\$25.30	
	1140C	1 L	\$39.70	
	1140E	1 gal	\$95.20	
Iodine, Verhoeff, Aqueous (Use: Verhoeff Elastic/Masson Trichrome procedure)	1209A	125 ml	\$16.00	
	1209B	250 ml	\$18.60	
Iodine, Weigert & Lugol, Aqueous (Use: Hucker-Twort Gram Stain, Verhoeff Elastic, McLetchie Trichrome and mercury pigment removal)	12092A	500 ml	\$31.60	
	12092B	1 L	\$44.50	
Lithium Carbonate, Saturated Aqueous (Use: Bluing Reagent)	12215B	1 L	\$36.80	
	12215C	1 gal	\$92.60	
Methenamine 3%, Aqueous (Use: Gomori Basement Membrane Stain (PAMM) and Urates)	12239A	250 ml	\$28.90	
	12239B	500 ml	\$50.70	
Mold Release Concentrate (Use: Help blocks release from embedding molds)	1246A	120 ml	\$14.00	H
	1246A	6 x 120 ml	\$50.40	H
Oxalic Acid 5%, Aqueous (Use: Bleach out melanin or phosphotungstic acid in the PTAH Stain and is stock for Gordon and Sweets Reticulum)	1293A	250 ml	\$39.40	
	1293B	500 ml	\$61.80	
Periodic Acid 0.5%, Aqueous (Use: PAS for glycogen or amoeba and in Gomori Basement Membrane Stain)	13308A	250 ml	\$31.80	
	13308B	500 ml	\$55.40	
Periodic Acid 1%, Aqueous (Use: Can be diluted to 0.5% for PAS or PAMM)	13304A	500 ml	\$59.50	
Phosphomolybdic - Phosphotungstic Acid, Aqueous (Use: Masson Trichrome Stain)	1332A	250 ml	\$72.30	
	1332B	500 ml	\$110.40	
	1332C	1 L	\$209.80	
Phosphomolybdic Acid 1%, Aqueous (Use: Klatskin Trichrome for liver)	13322A	500 ml	\$47.50	
	13322B	1 L	\$73.30	
Phosphomolybdic Acid 10%, Aqueous	13323A	500 ml	\$102.30	H
	13323B	1 L	\$177.60	H
Phosphotungstic - Acetic Acid, Aqueous (Use: With Gomori One-Step Trichrome)	13341A	500 ml	\$48.50	
	13341B	1 L	\$77.60	
Phosphotungstic Acid 2%, Alcoholic (Use: Trichrome, McLetchie)	13342B	500 ml	\$95.10	H
	13342C	1 L	\$135.90	H
Phosphotungstic Acid 4%, Alcoholic	13343B	500 ml	\$95.10	H
	13343C	1 L	\$135.90	H



PRODUCT NAME	PART	VOLUME	PRICE	
Phosphotungstic Acid 5%, Aqueous (Use: Trichrome, Masson Light Green)	13345A	250 ml	\$87.10	
	13345B	500 ml	\$129.50	
Picric Acid - Acetone 0.05% (Use: Brown and Hopps)	13351B	500 ml	\$64.20	H
	13351C	1 L	\$106.40	H
Picric Acid - Acetone 0.1% (Use: Gram Stains, e.g., Brown-Brenn)	1335A	500 ml	\$61.40	H
Picric Acid 7%, Alcoholic Stock (Use: Differentiator for Sulfated Alcian Blue stain and for formalin pigment removal)	1337A	500 ml	\$56.90	H
	1337B	1 L	\$67.60	H
	1337C	4 L	\$182.30	H
Picric Acid, Saturated Aqueous (Use: Bouin fixative)	1336A	500 ml	\$33.10	
	1336B	1 L	\$62.40	
	1336C	4 L	\$89.30	
Potassium Ferricyanide 1%, Aqueous (Use: Schmorl Melanin Stain)	13390A	250 ml	\$20.10	
	13390B	500 ml	\$36.10	
Potassium Ferricyanide 10%, Aqueous (Use: Schmorl Method for reducing substances)	13394A	500 ml	\$36.60	
Potassium Ferrocyanide 10%, Aqueous (Use: Gomori Mod. Iron Stain and stock for Colloidal Iron)	13392A	500 ml	\$48.10	
Potassium Permanganate 0.25%, Aqueous (Use: PTAH Stain and working solution for melanin bleaching)	133931A	250 ml	\$22.10	
	133931B	500 ml	\$30.00	
Potassium Permanganate 1%, Aqueous (Use: Gordon and Sweets Reticulum Stain and stock solution for melanin bleaching)	13393B	500 ml	\$38.20	
Propylene Glycol 100%, ACS (Use: Oil Red O in Propylene Glycol, Fat Stain)	13391B	250 ml	\$20.70	
	13391A	500 ml	\$38.90	
Propylene Glycol 85%, Aqueous (Use: Oil Red O in Propylene Glycol, Fat Stain)	133912A	250 ml	\$19.90	
	133912B	500 ml	\$25.20	
Rosin 10%, Alcoholic Stock (Use: Wolbach Giemsa)	13398A	500 ml	\$56.20	H
Sodium Bicarbonate 5%, Aqueous (Use: Original Brown and Brenn Gram Stain, and stock solution for Wolbach Giemsa)	13818A	500 ml	\$15.60	
Sodium Bisulfite 1%, Aqueous (Use: Grocott Methenamine Silver GMS Stain)	13821A	500 ml	\$23.70	
	13821B	1 L	\$37.70	
Sodium Borate 0.5%, Aqueous (Use: Rhodanine Stain for Copper)	13824A	250 ml	\$18.60	
	13824B	500 ml	\$23.10	
Sodium Borate 5%, Aqueous (Use: Standard Grocott Methenamine Silver, GMS procedure)	13826B	500 ml	\$33.30	

CHEMICALS & REAGENTS

PRODUCT NAME	PART	VOLUME	PRICE
Chemicals			
Sodium Metabisulfite 1%, Aqueous (Use: Working solution for GMS)	13853A	250 ml	\$16.10
	13853B	500 ml	\$22.60
	13853C	1 L	\$32.70
Sodium Sulfate 1%, Aqueous (Use: Sulfated Alcian Blue Stain)	1388A	500 ml	\$20.10
Sodium Thiosulfate 2%, Aqueous (Use: Working solution for GMS)	13888A	250 ml	\$22.70
	13888B	500 ml	\$26.80
Sodium Thiosulfate 2.5%, Aqueous (Use: Working solution for Gomori Basement Membrane Stain)	13889A	250 ml	\$14.80
	13889B	500 ml	\$17.50
Sodium Thiosulfate 3%, Aqueous (Use: Gomori Methenamine-Silver for Nitrates)	13887A	250 ml	\$15.80
	13887B	500 ml	\$26.00
Sodium Thiosulfate 5%, Aqueous (Use: Stock solution for Gomori Basement Membrane and GMS; Working solution for VVG Elastic, Gordon and Sweets and mercury pigment removal)	1389A	500 ml	\$28.90
	1389B	1 L	\$36.10
Sulfuric Acid 1%, Aqueous (Use: Fite Stain for Nocardia)	14012A	250 ml	\$14.10
	14012B	500 ml	\$20.10
	14012C	1 L	\$31.10
Uranyl Nitrate 1%, Aqueous (Use: Steiner Silver Impregnation Methods)	14036A	250 ml	\$85.40
	14036B	500 ml	\$164.90
	14036C	1 L	\$325.20
Water, Deionized ACS Reagent CLSI, Type 1	1408B	4 x 1 gal	\$24.30
	1408C	20 L cube	\$15.50



SLIDES & COVER SLIPS

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Non-Adhesive Slides				
Plain Approx. 144 slides/gross	75 x 25 x 1.0 mm	6210	gr	\$8.00
Frosted Approx. 144 slides/gross	75 x 25 x 1.0 mm One Side (Single)	6215	gr	\$8.75
Frosted Approx. 144 slides/gross	75 x 25 x 1.0 mm Both Sides (Double)	6216	gr	\$9.50
Colored End Approx. 144 slides/gross	White	6206W	gr	\$13.00
	Blue	6206B	gr	\$13.00
	Green	6206G	gr	\$13.00
	Pink	6206P	gr	\$13.00
	Yellow	6206Y	gr	\$13.00
Adhesive / Charged Slides				
APS				
APS Adhesive Microscope Slides Aminosilane coated glass slide with a hydrophobic surface. 75 x 25 x 1.0 mm. 1,000 slides/case	White	5078W	cs	\$265.00
	Blue	5078B	cs	\$265.00
	Green	5078G	cs	\$265.00
	Yellow	5078Y	cs	\$265.00
	Pink	5078P	cs	\$265.00
TruBond 360				
TruBond 360 Adhesive Microscope Slides Strong adhesion with a hydrophilic surface. 75 x 25 x 1.0 mm. 1,000 slides/case	White	5079W	cs	\$280.00
	Blue	5079B	cs	\$280.00
	Green	5079G	cs	\$280.00
	Yellow	5079Y	cs	\$280.00
	Pink	5079P	cs	\$280.00
TruBond 380				
TruBond 380 Adhesive Microscope Slides Very strong adhesion with a hydrophilic surface. 75 x 25 x 1.0 mm. 1,000 slides/case	White	5080W	cs	\$280.00
	Blue	5080B	cs	\$280.00
	Green	5080G	cs	\$280.00
	Yellow	5080Y	cs	\$280.00
	Pink	5080P	cs	\$280.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
TOMO				
TOMO Adhesive Microscope Slides Excellent adhesion with a hydrophilic surface. Designed to work on all autostainers. 75 x 25 x 1.0 mm. 1,000 slides/case	White	5081W	cs	\$280.00
	Blue	5081B	cs	\$280.00
	Green	5081G	cs	\$280.00
	Yellow	5081Y	cs	\$280.00
	Pink	5081P	cs	\$280.00
NewSilane				
NewSilane Very strong adhesion with a hydrophobic surface. 75 x 25 x 1.0 mm, approx. 72 slides/box	White	5070	1-10 bx	\$56.00/bx
		5070	11-19 bx	\$47.90/bx
		5070	20+ bx	\$41.90/bx
Diamond White				
Diamond White Charged Slides Charged surface. 75 x 25 x 1.0 mm, approx. 144 slides/gross, 10 gross/case	White	5085W	gr	\$40.20
		5085W	cs	\$321.00
	Aqua	5085A	cs	\$321.00
	Blue	5085B	cs	\$321.00
	Gold	5085D	cs	\$321.00
	Green	5085G	cs	\$321.00
	Lilac	5085L	cs	\$321.00
	Orange	5085O	cs	\$321.00
	Pink	5085P	cs	\$321.00
	Tan	5085T	cs	\$321.00
	Yellow	5085Y	cs	\$321.00
Swift Bond				
Swift Bond Strong aminosilane adhesion with a hydrophilic surface. 75 x 25 x 1.0 mm, 144 slides/gross, 10 gross/case	White	5075W	1-9 gr	\$37.00
			10+ gr	\$32.50
Superfrost Plus				
Superfrost Plus Charged surface. 75 x 25 x 1.0 mm, 10 gross/case	White	6203W	cs	\$722.60
Aminosilane 75 x 50 mm				
Aminosilane, 75 x 50 x 1.0 mm Aminosilane coated glass slide with a hydrophobic surface. 75 x 50 x 1.0 mm "Double Wide," approx. 144 slides/gross		5059F	gr	\$45.00

SLIDES & COVER SLIPS

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Premiere Aminosilane				
Premiere Charged Slides Aminosilane coated economy slide with a hydrophobic surface. 75 x 25 x 1.0 mm, 10 gross/case	White	6207W	cs	\$240.00
	Pink	6207P	cs	\$240.00
	Yellow	6207Y	cs	\$240.00
	Blue	6207B	cs	\$240.00
	Green	6207G	cs	\$240.00
Poly-L-Lysine				
Poly-L-Lysine, Newcomer Supply Brand Slides 75 x 25 x 1.0 mm, approx. 144 slides/gross	White	5010	1-9 gr	\$32.00
		5010	10-50 gr	\$30.00
		5010	51-99 gr	\$28.00
		5010	100+ gr	Call
	Blue	5013	1-9 gr	\$32.00
		5013	10-50 gr	\$30.00
		5013	51-99 gr	\$28.00
		5013	100+ gr	Call
Red Control Box Slides				
Red Control Box For Patient and Control tissue on the same slide. 75 x 25 x 1.0 mm, approx. 144 slides/gross	Red Box on back 1" Painted End PLUS	62220A	gr	\$93.20
	Red Box on back 3/4" Painted End PLUS	62221A	gr	\$91.10
	Red Box on back 1" Painted End	62222A	gr	\$90.50
	Red Box on front 3/4" Painted End PLUS	62223A	gr	\$97.00
	Red Box on front 3/4" Painted End	62224A	gr	\$87.10
Cover Slips & Mounting Mediums				
Cover Slips				
Premiere Brand Cover Slips Each case comes with 10 individual shrink sealed and desiccant 1 oz. boxes of the specified size cover slips.	18 mm x 18 mm #1	6230	10 oz/cs	\$25.00
	22 mm x 22 mm #1	62301	10 oz/cs	\$25.00
	22 mm x 22 mm #2	62302	10 oz/cs	\$25.00
	24 mm x 40 mm #1	62305	10 oz/cs	\$25.00
	24 mm x 50 mm #1	62306	10 oz/cs	\$25.00
	24 mm x 60 mm #1	62307	10 oz/cs	\$25.00
	55 mm x 50 mm #1	62308	50/pk	\$25.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE	
Mounting Mediums					
Permanent					
Permaslip Mounting Medium and Liquid Cover Slip. Fast drying, clear acrylic medium. Unique spout on 8 oz. cans controls drop by drop application.		6530B	1 to 11	\$35.00/ea	H
		6530B	12	\$279.60	H
Mount-Quick Xylene Soluble Mounting Medium. Fast drying and clear. UNIQUE AS A TISSUE TRANSFER MEDIUM. 30cc tube.		6270	each	\$17.50	
S Mounting Medium Compatible with xylene and aliphatic hydrocarbon xylene substitutes (e.g., SlideBrite). Dries clear and fast. Incompatible with film cover slips. 16 oz. bottle.		6750	each	\$46.00	H
		6750	12/cs	\$498.00	H
Limonene-Mount Permanent non-solvent based mounting medium made using limonene as a solvent instead of toluene or other petroleum pollutants.		51701A	100 ml	\$144.00	
		51701B	250 ml	\$295.00	
		51701C	6 x 100 ml	\$817.00	
O-Mount O-Mount is designed as a routine cover slipping mounting medium for organic insoluble chromogens, such as DAB. It also works well for routine H & E staining.		GBE02-18	18 ml	\$49.00	
Aqueous					
Advantage Permanent Mounting Medium preserves AEC, DAB, Fast Red and other dye chemistry stains.	Covers approx. 1,000 Slides	NB300	60 ml	\$295.00	
	Covers approx. 2,000 Slides	NB300A	125 ml	\$395.00	
Fluoro-Gel Aqueous based mounting medium designed for the preservation of fluorescence stained tissue sections.		56001A	20 ml	\$68.00	
		56001B	100 ml	\$210.00	
H&E Mount Unique, permanent mounting media designed for cover slipping H&E stains from water or alcohol.	Covers approx. 2,000 Slides	NB315	125 ml	\$325.00	
Mount-Quick Aqueous Aqueous Mounting Medium. Liquid at RT. Suitable for use with lipid sections, solvent soluble IHC chromogens and frozen sections. 30 cc, pH 7.		6271A	each	\$25.00	
Probe Mount Nuclear Fast Red stained specimens can be permanently mounted. Mount immunofluorescence and in-situ slides from water.	Covers approx. 1,000 Slides	NB316	60 ml	\$325.00	

SLIDES & COVER SLIPS

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Aqueous				
GB Mount Water-based mounting medium for frozen sections, immunoperoxidase and immunofluorescent procedures.		GBE01-18	18 ml	\$59.00
		GBE01-100	100 ml	\$292.00
Fluorescent Mounting Medium Aqueous mounting medium for preserving fluorescence of tissue and cell smears.		GBE18-18	18 ml	\$50.00
		GBE18-100	100 ml	\$279.00
Fluorescent Mounting Medium with DAPI (4,6-diamidino-2-phenylindole) Fluorescent mounting medium with DAPI is an aqueous mounting medium for preserving fluorescence of tissue and cell smears.		GBE19-18	18 ml	\$89.00
		GBE19-100	100 ml	\$427.00
Fluorescent Mounting Medium with PI (Propidium Iodide) Fluorescent mounting medium with PI is an aqueous mounting medium for preserving fluorescence of tissue and cell smears.		GBE20-18	18 ml	\$89.00
		GBE20-100	100 ml	\$427.00
Cover Glass Free				
Clear-Mount Permanent water-based mounting medium for Immunoperoxidase and immunoalkaline phosphatase IHC stains.		51921A	30 ml	\$72.00
		51921B	250 ml	\$275.00
ISH Mounting Medium Use for permanent preservation of tissue sections and cell smears with alkaline phosphatase chromogen, such as BCIP/NBT and red stain.		GBE21-18	18 ml	\$50.00
		GBE21-100	100 ml	\$279.00
Simpo-Mount Permanent water-based mounting medium for IHC stains. Simpo-Mount is a water-based mounting medium for immunohistology.		GBE03-18	18 ml	\$39.00
		GBE03-100	100 ml	\$130.00
Adhesive Solutions				
Aminoalkylsilane Slide Adhesive Shelf Life: 6 mos.		1007A	250 ml	\$28.30
		1007B	500 ml	\$37.40
		1007C	1L	\$50.30
Chrome Alum-Gelatin Adhesive (Use: Extra Sticking Power) For water bath. Store at 15-30°C.		1033A	250 ml	\$52.80
		1033B	500 ml	\$75.30
		1033C	1L	\$120.90
Egg Albumin Adhesive For direct slide application. Store at 4°C.		1058A	250 ml	\$41.60



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Haupt Gelatin Adhesive (Use: Gelatin and glycerin precombined for easy dissolving in water bath)		1151A	500 ml	\$45.50
Poly-L-Lysine Adhesive Stock Store at 4°C. Dilute 1 to 10.		1339A	100 ml	\$34.70

LAB ESSENTIALS & EQUIPMENT

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Laboratory Essentials				
Brushes				
Camel Hair Round Large Uses include: dusting microtome, cryostat and catching sections. Max brush diameter 3/4" x 1 1/2" long.		5107A	each	\$15.00
Camel Hair Round Small Uses include: preventing cryo sections from curling and laying out tissue in water bath. Brush diameter 3/16" x 13/16" long.		5107B	each	\$3.75
Camel Hair 1" Flat Uses include: general cleanup and removing air bubbles from water bath. Brush width 1" x 1" long.		5107C	each	\$7.50
Microtome & Cryostat Aligners				
Microtome Aligner - Universal Only one aligner needed for ALL models.		6257TO	each	\$250.00
Cryoaligner (Cryostat Alignment Tool) Will easily align both the tissue chuck and blade holder in perfect alignment within minutes. Can also be used as a heat extractor due to its stainless steel construction.		5190	each	\$199.00
Microtome Blades				
C.L. Sturkey Microtome Blades Select Teflon finish to create an edge that can produce superior sections with every cut. It is ideal for all general paraffin applications and cryostat use.	Low Profile 50 Blade	62514A	1-4	\$78.00
	Dispenser (Price per box)	62514A	5-9	\$76.00
		62514A	10+	\$74.00
	High Profile 50 Blade	62514D	1-4	\$80.00
	Dispenser (Price per box)	62514D	5-9	\$78.00
		62514D	10+	\$76.00
C.L. Sturkey Microtome Blades Extremus Coated with the hardest of the nitride coatings to suit the widest variety of cutting conditions. Extremely long lifespan.	Low Profile 50 Blade	62512A	1-4	\$82.00
	Dispenser Box (Price per box)	62512A	5-9	\$80.00
		62512A	10+	\$78.00
	Low Profile 200 Blade	62512C	each	\$304.00
	Dispenser			
	High Profile 50 Blade	62512D	1-4	\$82.00
	Dispenser (Price per box)	62512D	5-9	\$80.00
		62512D	10+	\$78.00
	High Profile 120 Blade	62512F	each	\$186.00
	Dispenser			



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE		
C.L. Sturkey Microtome Blades Gold Ultimate blade for cutting hard tissue and bone. Specialty titanium, nitride coated.	Low Profile 50 Blade Dispenser Box (Price per box)	62511A	1-4	\$74.00		
		62511A	5-9	\$72.00		
		62511A	10+	\$70.00		
	Low Profile 200 Blade Dispenser	62511C	each	\$264.00		
		High Profile 50 Blade Dispenser (Price per box)	62511D	1-4	\$74.00	
			62511D	5-9	\$72.00	
	62511D		10+	\$70.00		
	High Profile 120 Blade Dispenser	62511F	each	\$176.00		
		C.L. Sturkey Microtome Blades Silver Excellent all-purpose blade for paraffin and frozen tissue. Blades have a nitride and copolymer finish.	Low Profile 50 Blade Dispenser Box (Price per box)	62513A	1-4	\$78.00
				62513A	5-9	\$76.00
	62513A			10+	\$74.00	
	Low Profile 200 Blade Dispenser Box	62513C	each	\$278.00		
High Profile 50 Blade Dispenser (Price per box)		62513D	1-4	\$80.00		
		62513D	5-9	\$78.00		
	62531D	10+	\$76.00			
High Profile 120 Blade Dispenser	62513F	each	\$186.00			
	C.L. Sturkey Microtome Blades Monarch For general paraffin and cryostat use.	Low Profile 50 Blade Dispenser (Price per box)	62515A	1-4	\$94.00	
			62515A	5-9	\$92.00	
62515A			10+	\$90.00		
OCT						
OCT, PolarStat Frozen Tissue Medium Sections at -10°C and below. Freezes rapidly and dissolves without residue.	4 oz. dispenser bottle	6405	each	\$13.00		
		6405	12/cs	\$144.00		
Cryo Spray 60 Below Chills to -60°F. Contains no CFCs or HCFCs. Straw attached to can to localize spray onto blocks and fresh tissue.	10 oz. spray can	51912	each	\$16.80		
		51912	12/cs	\$168.00		
PAP Pens						
PAP Pen Liquid Blocker, Super Hydrophobic liquid creates a secure boundary on glass or plastic around the tissue specimen for IHC procedures. Fluid flows as a light blue color on the slide.		6505	each	\$63.00		
		6505	6/set	\$360.00		
PAP Pen Liquid Blocker, Mini Same as Super PAP pen in a smaller version with finer tip.		6506	each	\$46.00		
		6506	6/set	\$266.40		

LAB ESSENTIALS & EQUIPMENT

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
PAP Pens				
PAP Gnome Pen™ Liquid Blocker Hydrophobic colored liquid creates a secure boundary on glass slides around the tissue specimen for IHC procedures. It dries in less than one minute and can produce very thin boundaries.	Blue	6507B	each	\$74.00
	Green	6507G	each	\$74.00
	Orange	6507O	each	\$74.00
	Purple	6507P	each	\$74.00
	Red	6507R	each	\$74.00
	Yellow	6507Y	each	\$74.00
	All Color Pack (6 colors)	K6507X	6 pk	\$429.00
Parafilm M				
Parafilm M Sealing Film Self sealing film. Holds moisture loss to a minimum. Moldable and flexible.	50 mm x 75 m	6538A	each	\$41.00
	100 mm x 38 m	6538B	each	\$41.00
	100 mm x 75 m	6538C	each	\$69.00
Parafilm M Dispenser, ABS, Blue Store, dispense and cut parafilm. Made of ABS plastic. Holds one or two 2-inch rolls or one 4-inch roll.		6538DISP	each	\$49.00
Pens, Marker				
HistoTec Pen Ultra fine-tip, black ink pen that writes and dries quickly. Works on microscope slides and plastic cassettes. Withstands histo & cyto stains, xylene, alcohol and other histo chemistries.		5725	12/bx	\$27.00
Moist Mark Plus Pen Black, fine-tip ink pen designed for histo and cyto labs. Writes without smearing on both slides and plastic cassettes.		6260	10/pk	\$31.00
Pipettes				
Micropipettes & Tips				
Micropipette: Adjustable Volume/Tip Ejector Super-soft elastomer grip that minimizes fatigue of hands and fingers. Comes with a certificate of calibration and a recalibration tool kit.	0.5 - 10 µl	6555A	each	\$177.00
	2 - 20 µl	6555B	each	\$177.00
	5 - 50 µl	6555C	each	\$177.00
	10 - 100 µl	6555D	each	\$177.00
	20 - 200 µl	6555E	each	\$177.00
	100 - 1000 µl	6555F	each	\$177.00
	1000 - 5000 µl	6555G	each	\$177.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Pipette Tips for Diamond Pro Pipettors	0.1 - 20 µl	6556	1,000/bag	\$18.00
		6556R	96/rk, 10 rk/bx	\$47.00
	1 - 200 µl	6557	1,000/bag	\$20.00
		6557R	200/rk, 5 rk/bx	\$37.00
	100 - 1000 µl	6558	1,000/bag	\$24.00
		6558R	200/rk, 5 rk/bx	\$40.00
	1000 - 5000 µl	6559	250/bag	\$39.00
	Transfer Pipettes			
Pipets, Transfer		6550	500/bx	\$18.00
Polyethylene pipet with large bulb. Pipet graduated from 0.5 to 3.0 ml.				
Slide Boxes				
Microscope Slide Box - Small	White	62421A	each	\$2.50
	Red	62421B	each	\$2.50
	Yellow	62421C	each	\$2.50
	Blue	62421D	each	\$2.50
	Green	62421E	each	\$2.50
Microscope Slide Box - Medium	White	62422A	each	\$5.80
	Red	62422B	each	\$5.80
	Yellow	62422C	each	\$5.80
	Blue	62422D	each	\$5.80
	Green	62422E	each	\$5.80
Microscope Slide Box - Large	White	6243A	each	\$7.00
	Red	6243B	each	\$7.00
	Yellow	6243C	each	\$7.00
	Blue	6243D	each	\$7.00
	Green	6243E	each	\$7.00
	Gray	6243F	each	\$7.00
	Black	6243G	each	\$7.00
Microscope Slide Box - Tall (75 mm x 50 mm slides)	1 to 9	6244	each	\$18.00
	10 or more	6244	each	\$16.00
Identification card on inside lid. Stackable with other boxes. Holds up to 100 slides.				

LAB ESSENTIALS & EQUIPMENT

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Slide Draining Trays				
Slide Draining Tray, Polyethylene Holds slides in a vertical position for drying.	15 Slide Capacity	6842B	each	\$30.00
	40 Slide Capacity	6842C	each	\$33.00
	78 Slide Capacity	6842D	each	\$41.00
Slide Mailers				
Slide Mailers, Cardboard Made of sturdy cardboard to protect the slides. Thumb groove for easy slide removal.	1 Place Mailer	68412	100/pk	\$17.00
		68412	2,000/cs	\$225.00
	2 Place Mailer	68413	50/pk	\$14.00
		68413	1,000/cs	\$195.00
	4 Place Mailer	68414	25/pk	\$13.00
		68414	200/cs	\$76.00
Slide Mailer, Plastic, 1 Place Polyethylene mailers have a sturdy resealable clasp. Ideal for shipping specimen slides.		68411	50/pk	\$11.50
		68411	2,000/cs	\$310.00
Slide Mailer, Plastic, 5 Place Polypropylene, rectangular slide box with hinged flip top lid. 5 slide capacity.		6841	100/bx	\$24.00
		6841	800/cs	\$160.00
Slide Mailer, Plastic, 5 Place, Assorted Colors (Five of each color: Red, Yellow, Blue, Purple, and Orange)		6841C	25/pk	\$23.00
Slide Presentation Folders				
Slide Folder, Cardboard Cardboard slide trays hold standard (75 x 25 x 1.0 mm) slides. The horizontal compartments are recessed protecting the specimens when the cover is closed.	10 Place	6839	1-24	\$5.50
		6839	25-99	\$4.60
		6839	100+	\$3.90
	20 Place	6840	1-24	\$6.00
		6840	25-99	\$5.00
		6840	100+	\$4.25
	30 Place	68401	1-24	\$12.00
		68401	25-99	\$9.50
		68401	100+	\$8.00
Slide Folder Rack for 20 Place Slide Folder Store multiple 20 place slide trays for easy retrieval. Stackable.		6840RK	each	\$31.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Slide Folder, Plastic Opaque, plastic back with hinged covers that will fold under folder. Clear plastic cover. Holds 20 standard slides.	Blue	68402B	each	\$10.00
		68402B	10/pk	\$86.00
	Green	68402G	each	\$10.00
		68402G	10/pk	\$86.00
	Gray	68402GY	each	\$10.00
		68402GY	10/pk	\$86.00
	White	68402W	each	\$10.00
		68402W	10/pk	\$86.00
	Yellow	68402Y	each	\$10.00
		68402Y	10/pk	\$86.00
Slide Tray, Plastic Stackable trays so that slide will not touch tray above. Holds 20 standard slides.	White	6845W	each	\$6.00
Slide Staining Trays - IHC				
Slide Moisture Chamber 10 slide chambers for maintaining moisture and watching the reaction. Stackable, heavy duty plastic construction.	Clear	68431A	each	\$49.20
	Black	68432A	each	\$49.20
SlideShow™ 30	Black Lid, holds up to 30 slides	6844-30BL	each	\$235.00
	Clear Lid, holds up to 30 slides	6844-30CL	each	\$235.00
SlideShow™ 20	Clear Lid, holds up to 20 slides	6844-20CL	each	\$170.00
SlideShow™ 10	Clear Lid, holds up to 10 slides	6844-10CL	each	\$155.00
StainTray 20 Place	Black Lid, holds up to 20 slides	6847-20BL	each	\$182.00
	Clear Lid, holds up to 20 slides	6847-20CL	each	\$182.00
	Clear Lid Only	6847-20LIDC	each	\$67.00
	Black Lid Only	6847-20LIDB	each	\$67.00
StainTray 10 Place	Black Lid, holds up to 10 slides	6848-10BL	each	\$156.00
	Clear Lid, holds up to 10 slides	6848-10CL	each	\$156.00
	Clear Lid Only	6848-10LIDC	each	\$47.00
	Black Lid Only	6848-10LIDB	each	\$47.00
Slide Staining Trays - Disposable Tray holds up to 8 slides.		6846	4 bases & 1 lid	\$39.00

LAB ESSENTIALS & EQUIPMENT

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Specimen Chucks				
Specimen Chuck With black o-rings.	Round, 25 mm	6870	each	\$45.00
	Round, 40 mm	6872	each	\$68.00
	Extender, 1 5/8" x 1 1/8"	6871	each	\$60.00
Specimen Chuck Colored Set, 25 mm w/Orientation Mark Set includes: 8-25 mm specimen chucks, 1 of each color (Pink, Yellow, Green, White, Black, Purple, Blue, and Red)		6866S	set	\$535.00
Specimen Chuck w/Colored O-ring and Orientation Mark, Round 25 mm	Red	6867R	each	\$40.00
	White	6867W	each	\$40.00
	Green	6867G	each	\$40.00
	Yellow	6867Y	each	\$40.00
	Pink	6867P	each	\$40.00
	Blue	6867B	each	\$40.00
Specimen Chuck w/Colored O-ring and Orientation Mark, Round 40 mm	Red	6868R	each	\$75.00
	Blue	6868B	each	\$75.00
	Pink	6868P	each	\$75.00
	White	6868W	each	\$75.00
	Green	6868G	each	\$75.00
	Yellow	6868Y	each	\$75.00
Specimen Chuck w/Colored O-ring and Orientation Mark, Kit Includes: 10 specimen chucks (6-30 mm, 2-25 mm, 2-40 mm), teasing needle, and angular brush		6873S1	kit	\$365.00
Specimen Chuck w/Colored O-ring and Orientation Mark, Round 55 mm	White	6869W	each	\$78.00
	Red	6869R	each	\$78.00
	Dark Purple	6869U	each	\$78.00
	Green	6869G	each	\$78.00
	Blue	6869B	each	\$78.00
	Yellow	6869Y	each	\$78.00

Staining Dishes & Racks

Coplin Jars

Coplin Jar, Plastic, White Can hold up to 10 slides and is microwave safe. Pricing is per unit.	5184	1-5	\$8.00
	5184	6-11	\$7.00
	5184	12+	\$5.00
Coplin Jar, Glass 5 horizontal slots. Includes glass cover.	5185	each	\$12.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Glass Staining Dishes & Racks				
Staining Dish, Glass 10 horizontal slots hold 10 to 18 standard slides. Includes glass cover. Not for use with slide racks.		6815	each	\$14.00
Staining Dishes The dishes are made of a heavy glass with snug fit lids. (Designed to work with part #6817 racks)	For 20 slide rack	6816A	each	\$20.40
	For 30 slide rack	6816B	each	\$22.80
	For 60 slide rack	6816C	each	\$36.00
Staining Racks Stainless steel racks with handles. (Designed to work with part # 6816 dishes)	Holds 20 slides horizontally	6817A	each	\$16.80
	Holds 30 slides vertically	6817B	each	\$28.00
	Holds 60 slides vertically	6817C	each	\$37.20
Easy Dip™ Staining System				
EasyDip™ Slide Staining Jar Resistant to most chemicals and dyes. Withstands temperatures up to 121°C. Jars all have attached hinged lids. (Jars are not resistant to: phenol, iodine, ferric chloride, hydrochloric acids or hydrogen peroxide. Rack sold separately)	Blue	5300B	each	\$11.00
		5300B	6/cs	\$58.20
	Green	5300G	each	\$11.00
		5300G	6/cs	\$58.20
	Pink	5300P	each	\$11.00
		5300P	6/cs	\$58.20
	White	5300W	each	\$11.00
		5300W	6/cs	\$58.20
	Yellow	5300Y	each	\$11.00
		5300Y	6/cs	\$58.20
EasyDip™ Slide Staining Rack Holds up to 12 slides vertically. Writing area of slide will not be stained. Rack cannot be microwaved.	Gray	5300RK	each	\$12.00
		5300RK	6/cs	\$63.60
EasyDip™ Slide Staining Kit Contains: 5 jars (one of each color) and 1 rack.		K5300	each	\$60.00
EasyDip™ Anodized Aluminum Jar Rack Holder Securely holds up to six Easy Dip jars and prevents them from tipping.		5300JRK	each	\$71.00

LAB ESSENTIALS & EQUIPMENT

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Sterile Instrument Packs				
Instrument Pack, Sterile Each envelope contains a stainless steel 4" forcep, a plastic handled scalpel with a #11 blade and a pointed scissors.		5810	12/pk	\$101.00
Storage Systems for Slides and Blocks				
Histology Slide & Tissue Block Storage System Cardboard construction and easy assembly. (Block Storage Box has 8 drawers each, Slide Storage Box has 4 drawers each)	Slide Storage System	6875	10 bx/cs	\$72.00
	Block Storage System	6877	10 bx/cs	\$87.00
Blocks Out/In Cards Mark when a block has been removed for easy record keeping until it is returned to the file system.		6878	1,000/bx	\$33.00
Slides Out/In Cards Mark when a slide has been removed for easy record keeping until it is returned to the file system.		6876	1,000/bx	\$33.00
Index Tab Cards Allow for easy categorizing and partitioning of slides in the storage file drawers.		6879	500/bx	\$27.00
Storage File From/To Stickers Used to place on the recording area of the front panel to reuse the box. The sticker can also be used on the front of each drawer.		6880	500/roll	\$27.00
Timers				
Pocket Timer 4 Channel w/Clock Four independent timing channels with visual and audible alarms. Includes: spring clip, magnetic back, and batteries.		6886	each	\$25.00
Digital Timer, Sand Count Up/Down Audible and visual alarms. One side displays count up/down timer, the other side displays the clock in 12- or 24-hour Single Channel.		6887	each	\$27.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Xylene Pads				
Hydrosorb-X Water absorbent packets designed to fit in the bottom of your clearing agent dish (manual or auto). Extends life of clearant and reduces disposal costs.	3 1/4" x 1 1/2"	STHXP-325-150	10/bx	\$70.00
	3 1/4" x 2 1/2"	STHXP-325-250	10/bx	\$70.00
	3 1/2" x 1 1/4"	STHXP-350-125	10/bx	\$70.00
	3 1/2" x 1 1/2"	STHXP-350-150	10/bx	\$70.00
	3 1/2" x 2 3/4"	STHXP-350-275	10/bx	\$70.00
	4" x 1 1/2"	STHXP-400-150	10/bx	\$70.00
	4" x 2"	STHXP-400-200	10/bx	\$70.00
	4" x 2 1/2"	STHXP-400-250	10/bx	\$70.00
	4 1/2" x 1 1/2"	STHXP-450-150	10/bx	\$70.00
	5" x 1 1/2"	STHXP-500-150	10/bx	\$70.00
	3 3/4" x 2 3/4"	STHXP-375-275	10/bx	\$75.00
	3 3/4" x 3"	STHXP-375-300	10/bx	\$75.00
	3" x 4"	STHXP-300-400	10/bx	\$75.00
	4 1/4" x 3"	STHXP-425-300	10/bx	\$75.00
	5" x 3"	STHXP-500-300	10/bx	\$75.00
	4" x 5"	STHXP-400-500	10/bx	\$80.00
	7" x 1 1/2"	STHXP-700-150	10/bx	\$75.00

Laboratory Equipment

Alcohol Recycling

Gravity Recycling Alcohol Cartridge (GRAC)

For Histology and Cytology labs that are only recycling alcohol from a staining setup. Removes stains, dyes, and small particulates from reagent alcohol.

STGRAC-100 each \$425.00 H

Gravity Recycling Alcohol Cartridge, Tissue Processor (GRAC-TP)

Ideal for Histology labs that are recycling from both staining setups and tissue processors. Removes stains, dyes, lipids, salts, and proteins from reagent alcohol.

STGRAC-TP each \$565.00 H

LAB ESSENTIALS & EQUIPMENT

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Alcohol Recycling				
One Gallon Collection System For small to medium labs (up to 40,000 slides per year). Includes: 1-one gallon "Dirty" container, 1-one gallon "Recycled" container, 1 filter sock, 1 drip tray, and instruction manual.		STCS-1	set	\$105.00
Two Gallon Collection System For medium to large labs (40,000+ slides per year). Includes: 1-two gallon "Dirty" container, 1-two gallon "Recycled" container, 1 filter sock, 1 drip tray, and instruction manual.		STCS-2	set	\$345.00
Five Gallon Collection System For labs recycling both processor and tissue processing alcohols. Includes: 1-five gallon "Dirty" container, 1-five gallon "Recycled" container, 1 filter sock, 1 drip tray, and instruction manual.		STCS-5	set	\$435.00
Five Gallon Additional container for "Clean" alcohol storage.		STCS-5RC	each	\$260.00
Hydrometer Kit Monitor the concentration of the recycled alcohol.		STHK-250	kit	\$65.00
Filter Sock, 14", 1 micron Replacement pre-filter for removal of cell debris and tissue fragments.		STFS-14A	each	\$23.50
Centrifuges				
Mini Centrifuge XC-4610K Perfect for quick spinning of small samples. It comes standard with 2 angle rotors and various speed options.		5134	each	\$250.00
Low Speed Laboratory Centrifuge XC-L4 Main unit does not come with any rotors. Must order at least one of the four rotor options.	Low Speed XC-L4 Centrifuge	5135A	each	\$610.00
	Fixed Angle Rotor for XC-L4 10 ml X 18 (4000 r/min), RCF: 2220 X g	5135B	each	\$230.00
	Fixed Angle Rotor for XC-L4 10 ml X 24 (4000 r/min), RCF: 2220 X g	5135C	each	\$230.00
	Fixed Angle Rotor for XC-L4 20 ml X 12 (4000 r/min), RCF: 2220 X g	5135D	each	\$230.00
	Swing Rotor for XC-L4 10 ml X 8 (4000 r/min), RCF: 2220 X g	5135E	each	\$515.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Low Speed Laboratory Centrifuge XC-L5A Main unit does not come with any swing rotors. Must order at least one of the six swing rotor options.	XC-L5A Model Centrifuge	5133A	each	\$2,212.00
	Swing Rotor for XC-L5A 7 ml X 48 (5000 r/min), RCF: 5030 X g	5133B	each	\$713.00
	Swing Rotor for XC-L5A 250 ml X 4 (5000 r/min), RCF: 5030 X g	5133C	each	\$993.00
	Swing Rotor for XC-L5A 50 ml X 8 (5000 r/min), 5030 X g	5133D	each	\$511.00
	Swing Rotor for XC-L5A 10 ml X 32 (5000 r/min), RCF: 5030 X g	5133E	each	\$568.00
	Swing Rotor for XC-L5A 500 ml X 4 (4000 r/min), RCF: 3220 X g	5133F	each	\$1,362.00
	Swing Rotor for XC-L5A 8 X 96 Hole (3000 r/min) RCF: 2120 X g	5133G	each	\$1,079.00
High Speed Microhematocrit Centrifuge XC-H12 Capillary Rotor included.		5136	each	\$1,540.00
Cryo Bath				
Cryo Bath for Histology Working temperature: fixed at -60°C.		5191	each	\$3,169.00
Chuck Holder for Cryo Bath For additional capacity. (One chuck holder is included with the Cryo Bath unit.)		5191CH	each	\$45.00
Cryostat Freeze Plate				
Cryostat Freeze Plate Working temperature of surface - 30°F.		51911	each	\$1,920.00
Keyboards & Mouse - Disinfectable				
Keyboard Water-resistant and disinfectable.		5843	each	\$199.00
Keyboard with Touchpad Mouse Water-resistant and disinfectable. Touch pad mouse built into the keyboard.		5844	each	\$259.00
Mouse Water-resistant and disinfectable.		5845	each	\$96.00

LAB ESSENTIALS & EQUIPMENT

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Microtomes				
Microtome, Manual TN5000	Manual	62531	each	\$8,350.00
Microtome, Semi-Automatic TN5500	Semi-Automatic	62532	each	\$8,500.00
Paraffin Dispenser				
Paraffin Dispenser Capacity 2.3 gallon.		6527	each	\$799.00
Paraffin Wax Trimmers				
Paraffin Block Trimmer, Programmable, TN-100 Programmable to set desired temperature at startup.		6529	each	\$452.00
Paraffin Wax Trimmer Two temperature settings: Low - approx. 80°F and High - approx. 90°F.		6528	each	\$365.00
Paraffin Wax Trimmer Disposable Liners for Tray A plastic tray liner for easy clean-up and disposal of the paraffin. (To be used with part #6528)		6528L	50/pk	\$23.50
Slide Warmers				
Slide Warmer To dry paraffin tissue sections on slides.	10" x 7" surface (Capacity - 23 slides)	6258	each	\$235.00
	25" x 8" surface (Capacity - 66 slides)	6259	each	\$300.00
Slide Warmer with Cover To dry paraffin tissue sections on slides. Heated surface can be covered with lid.	10" x 7" surface (Capacity - 23 slides)	6262	each	\$340.00
	25" x 8" surface (Capacity - 66 slides)	6261	each	\$400.00
Slide Warmer Step Up To dry paraffin tissue sections on slides. Can accommodate twice as many slides as a flat unit of the same size.		6259STEP	each	\$415.00
Stirring Hot Plate				
Stirring Hot Plate Has the ability to heat and stir simultaneously or perform each function separately.		6825	each	\$285.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Tissue Embedding Center				
Tissue Embedding Center, Programmable TN-1500 Used to embed pathological tissue and to cool tissue after embedding.		5310	each	\$8,950.00
Water Baths - General Purpose				
2 Liter Digital Laboratory Water Bath Digital control and programmability.		6955A	each	\$335.00
5 Liter Digital Laboratory Water Bath Digital control and programmability.		6955C	each	\$404.00
10 Liter Digital Laboratory Water Bath Digital control and programmability.		6955D	each	\$482.00
20 Liter Digital Laboratory Water Bath Digital control and programmability.		6955E	each	\$580.00
28 Liter Digital Laboratory Water Bath Digital control and programmability.		6955G	each	\$609.00
Water Baths - Tissue Flotation				
Tissue Flotation Bath - Digital Low profile design with a removable glass basin and is temperature controlled. The basin is illuminated from the side.		6950	each	\$600.00
Tissue Flotation Bath - Digital with Large Lighted Basin TN-450 Large bath with digital control and removable glass basin. LED light illuminates both front and back of bowl.		6952	each	\$1,325.00
Tissue Flotation Bath - Economy Analog controlled temperature.		6951	each	\$199.00

GROSSING & EMBEDDING

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Base Molds & Tampers				
HistoPress Tissue Embedding Tampers Tool for flattening tissue sections in an embedding mold.	1/2" x 1/2"	5103A	2/pk	\$59.00
	1/2" Diameter	5103B	2/pk	\$59.00
	3/4" x 9/16"	5103C	2/pk	\$59.00
	3/4" x 1 1/8"	5103D	2/pk	\$69.00
Base Molds, Stainless Steel Reusable stainless steel base molds.	7 x 7 x 5 mm	5104A	12/pk	\$34.00
	15 x 15 x 5 mm	5104B	12/pk	\$34.00
	24 x 24 x 5 mm	5104C	12/pk	\$34.00
	30 x 24 x 5 mm	5104D	12/pk	\$34.00
	37 x 24 x 5 mm	5104E	12/pk	\$34.00
Base Mold, Jumbosette, Stainless Steel Reusable stainless steel base molds.	66 x 45 x 15 mm	5104F	1-3 each	\$28.00
		5104F	4-6 each	\$25.50
		5104F	7-9 each	\$23.00
		5104F	10+ each	\$20.50
Base Molds, Disposable Designed for easy block release.	7 x 7 x 5 mm	5105A	500/cs	\$63.00
	15 x 15 x 5 mm	5105B	500/cs	\$63.00
	24 x 24 x 5 mm	5105C	500/cs	\$63.00
	30 x 24 x 5 mm	5105D	500/cs	\$63.00
	37 x 24 x 5 mm	5105E	500/cs	\$63.00
	37 x 24 x 10 mm	5105F	500/cs	\$63.00

Cassettes

Premiere Cassettes

Cassettes, Processing and Embedding Covers detached and packaged separately in box.	Blue	5115B	500/bx	\$42.50
	Gray	5115GY	500/bx	\$42.50
	Green	5115G	500/bx	\$42.50
	Lilac	5115L	500/bx	\$42.50
	Orange	5115O	500/bx	\$42.50
	Pink	5115P	500/bx	\$42.50
	Tan	5115T	500/bx	\$42.50
	White	5115W	500/bx	\$42.50
	Yellow	5115Y	500/bx	\$42.50
Cassettes, Biopsy Covers detached and packaged separately in box.	Blue	5116B	500/bx	\$42.50
	Green	5116G	500/bx	\$42.50
	White	5116W	500/bx	\$42.50
	Yellow	5116Y	500/bx	\$42.50



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Simport Cassettes				
Cassettes, Unisette (Processing & Embedding) Snap-latch and hinge-lock lid to the base.	Aqua	5124A	500/bx	\$44.00
	Blue	5124B	500/bx	\$44.00
	Green	5124G	500/bx	\$44.00
	Lilac	5124L	500/bx	\$44.00
	Pink	5124P	500/bx	\$44.00
	White	5124W	500/bx	\$44.00
	Yellow	5124Y	500/bx	\$44.00
Cassettes, Unisette Biopsy (Processing & Embedding) Snap-latch and hinge-lock lid to the base.	Aqua	5125A	500/bx	\$44.00
	Blue	5125B	500/bx	\$44.00
	Green	5125G	500/bx	\$44.00
	Lilac	5125L	500/bx	\$44.00
	Pink	5125P	500/bx	\$44.00
	White	5125W	500/bx	\$44.00
	Yellow	5125Y	500/bx	\$44.00
Cassettes, Slimsette (Processing & Embedding) Recessed lid, one-piece cassette, base and lid are attached via hinge.	Aqua	5126A	500/bx	\$65.00
	Blue	5126B	500/bx	\$65.00
	Green	5126G	500/bx	\$65.00
	Lilac	5126L	500/bx	\$65.00
	White	5126W	500/bx	\$65.00
	Yellow	5126Y	500/bx	\$65.00
Cassettes, Slimsette Biopsy (Processing & Embedding) Recessed lid, one-piece cassette, base and lid are attached via hinge.	Aqua	5127A	500/bx	\$65.00
	Blue	5127B	500/bx	\$65.00
	Green	5127G	500/bx	\$65.00
	Lilac	5127L	500/bx	\$65.00
	White	5127W	500/bx	\$65.00
	Yellow	5127Y	500/bx	\$65.00
Cassettes, Micromesh Biopsy No compartments. Lid Attached.	White	5120W	250/bx	\$59.00
	Yellow	5120Y	250/bx	\$59.00
	Pink	5120P	250/bx	\$59.00
	Lilac	5120L	250/bx	\$59.00
	Blue	5120B	250/bx	\$59.00
	Tan	5120T	250/bx	\$59.00
	Orange	5120O	250/bx	\$59.00
	Green	5120G	250/bx	\$59.00
Cassettes, Macrosette Specially designed for larger tissue specimens (10 mm deep). Hinged lid attached. Dimensions: 40.1 x 28.5 x 13 mm	White	5118W	250/bx	\$63.00
Cassettes, Jumbosette (Processing & Embedding) Dimensions: 75 x 52 x 18 mm Sold by each.	White	5128	1-9	\$4.00
		5128	10-24	\$3.60
		5128	25+	\$3.30

GROSSING & EMBEDDING

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Simport Cassettes				
Cassettes, Micromesh Biopsy (4 Compartments) Four square compartments each measuring 13 x 13 mm. Lid attached.	Blue	5119B	250/bx	\$59.00
	Green	5119G	250/bx	\$59.00
	Pink	5119P	250/bx	\$59.00
	White	5119W	250/bx	\$59.00
	Lilac	5119L	250/bx	\$59.00
	Yellow	5119Y	250/bx	\$59.00
Cassettes, Microsette Biopsy (6 Compartments) Six compartments each measuring 7 x 12 mm. Lid attached.	Blue	5117B	250/bx	\$59.00
	Green	5117G	250/bx	\$59.00
	Pink	5117P	250/bx	\$59.00
	White	5117W	250/bx	\$59.00
	Yellow	5117Y	250/bx	\$59.00
Cassettes, Slimsette in QuickLoad Sleeves (Processing & Embedding) 75 cassettes/sleeve. 10 sleeves/case.	White	5121W	750/cs	\$213.00
	Aqua	5121A	750/cs	\$213.00
	Pink	5121P	750/cs	\$213.00
	Green	5121G	750/cs	\$213.00
	Yellow	5121Y	750/cs	\$213.00
	Blue	5121B	750/cs	\$213.00
	Tan	5121T	750/cs	\$213.00
	Gray	5121GY	750/cs	\$213.00
	Peach	5121PE	750/cs	\$213.00
	Lilac	5121L	750/cs	\$213.00
	Orange	5121O	750/cs	\$213.00
Cassettes, Slimsette in QuickLoad Sleeves Biopsy 75 cassettes/sleeve. 10 sleeves/case.	White	5122W	750/cs	\$213.00
	Aqua	5122A	750/cs	\$213.00
	Pink	5122P	750/cs	\$213.00
	Green	5122G	750/cs	\$213.00
	Yellow	5122Y	750/cs	\$213.00
	Blue	5122B	750/cs	\$213.00
	Tan	5122T	750/cs	\$213.00
	Gray	5122GY	750/cs	\$213.00
Cassettes, Micromesh in QuickLoad Sleeves Biopsy (Processing & Embedding) 75 cassettes/sleeve. 10 sleeves/case.	White	5123W	750/cs	\$236.00
	Aqua	5123A	750/cs	\$236.00
	Pink	5123P	750/cs	\$236.00
	Green	5123G	750/cs	\$236.00
	Yellow	5123Y	750/cs	\$236.00
	Blue	5123B	750/cs	\$236.00
	Tan	5123T	750/cs	\$236.00
Biopsy Sponge Pads & Wraps				
Biopsy Sponge Pads for Cassettes Holds biopsies in place. Dimensions: 30.2 x 25.4 x 2 mm 10,000/case.		5106A	10,000/cs	\$223.00



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Biopsy Sponge Pads for Slimsette Cassettes Holds biopsies in place. Dimensions: 27.4 x 25.4 x 2 mm 10,000/case.		5106B	10,000/cs	\$223.00
Biopsy Wraps, Biosavers Biopsy wrap paper. Size: 2" x 3"	White	5550	1,000 sheets/bx	\$40.00
	Pink	5551	1,000 sheets/bx	\$40.00
	Blue	5552	500 sheets/bx	\$22.00
Baskets for Tissue Cassettes				
Cassette Basket, Small Holds up to 50 cassettes. Dimensions: 5 1/2" x 4 1/2" x 1 1/2"		5112A	each	\$270.00
Cassette Basket, Medium Holds up to 100 cassettes. Dimensions: 9" x 5 1/2" x 1 1/2"		5112B	each	\$345.00
Cassette Basket, Deluxe Holds up to 150 cassettes. Dimensions: 8 3/4" x 6" x 2"		5112C	each	\$560.00
Grossing Tools				
Dissecting Mat				
Dissecting Mat, Flexible, Polypropylene Metric grids to easily and quickly make measurements.	8" x 10"	5218A	10/pk	\$18.20
	12" x 18"	5218B	10/pk	\$39.10
Forceps				
Standard & Unique Forceps				
Standard Forceps, 11.5 cm		5511	each	\$9.90
Standard Forceps, 16 cm		5512	each	\$13.90
Standard Forceps, 20 cm		5513	each	\$15.90
Standard Embedding, 16 cm, Straight, Standard Grade, Smooth Tip Forceps		5538	each	\$17.80
Custom Embedding, 18 cm, Straight, Standard Grade Forceps		5540	each	\$15.20

GROSSING & EMBEDDING

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Standard & Unique Forceps				
Custom Embedding, 13 cm, Straight, Standard Grade Forceps		5535	each	\$13.90
Custom Embedding, 13 cm, Angled Curve, Standard Grade Forceps		5535B	each	\$15.30
Custom Embedding, 13 cm, Curved, Standard Grade Forceps		5536	each	\$15.30
Custom Embedding, 15 cm, Angled Curve, Standard Grade Forceps		5536B	each	\$15.80
Custom Embedding, 15 cm, Curved, Standard Grade Forceps		5536C	each	\$15.80
Custom Embedding, 15 cm, Straight, Standard Grade Forceps		5536D	each	\$13.20
Custom Embedding, 18 cm, Curved, Standard Grade Forceps		5537	each	\$16.60
Dissecting, Delicate, 12 cm, Curved, Standard Grade Forceps		5539	each	\$13.90
Ergonomic, 14 cm, 1 mm Serrated Tip, Premium Grade Forceps		5541A	each	\$53.00
Ergonomic, 14 cm, 1 mm Smooth Tip, Premium Grade Forceps		5541B	each	\$53.00
Ergonomic, 18 cm, 1:2 Teeth, Straight, Premium Grade Forceps		5542A	each	\$55.00
Ergonomic, 18 cm, 1 mm Serrated Tip, Premium Grade Forceps		5542B	each	\$55.00
Ergonomic, 18 cm, 1 mm Smooth Tip, Premium Grade Forceps		5542C	each	\$55.00
Ergonomic, 18 cm, 1 mm Smooth Tip, Angled, Premium Grade Forceps		5542D	each	\$55.00
Russian, 15 cm, Standard Grade Forceps		5543	each	\$23.80
12 cm Fine Plastic Tipped Forceps		5502	10/pk	\$18.50
15 cm Bent Fine Plastic Tipped Forceps		5503	10/pk	\$18.50
Flexible Plastic Grips Forceps Hold glass slides with ease.		5500	10/pk	\$22.00
Sterile Forceps Individually wrapped, stainless steel, 4 1/2".		5509A	25/pk	\$60.00
		5509B	100/pk	\$207.50



PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Lighted Forceps				
Each kit includes: one instrument, power pack and battery charger.				
13 cm, Smooth, Fine-Tip Embedding Forceps		5515	kit	\$250.00
15 cm, Smooth, Fine-Tip Embedding Forceps		5516	kit	\$250.00
18 cm, Smooth, Fine-Tip Embedding Forceps		5517	kit	\$250.00
14 cm, Smooth, Fine-Tip, Ergonomic Embedding Forceps		5518	kit	\$270.00
14 cm, Smooth, Fine Tip, Ergonomic Embedding Forceps (Tool Only)		5518IO	each	\$145.00
17.8 cm, Smooth, Fine-Tip, Ergonomic Embedding Forceps		5521	kit	\$275.00
16 cm, Blunt, Serrated Grossing Forceps		5519	kit	\$250.00
20 cm, Blunt, Serrated Grossing Forceps		5520	kit	\$250.00
Forceps Jars				
Forceps Jar - Stainless Steel Weighted bottom. Dimensions: 2" opening, 4 1/2" base diameter, 5" tall		5514B	each	\$58.00
Forceps Jar - Heated Three well heater. Temperature range of 180°F to 200°F.		5514	each	\$550.00
Scalpels				
Scalpels, Flat Edge, Standard Grade, 3.7 cm		5545	each	\$15.90
Scalpels, Finger Scalpels	Blade #10	6810A	100/cs	\$103.60
	Blade #22	6810B	100/cs	\$113.00
	Blade #60	6810C	100/cs	\$122.80
Scalpels, Stealth Sterile Safety Disposable and sterile, one-handed easy "slide & click" design.	Blade #22	6820A	100/cs	\$126.80
	Blade #60	6820B	100/cs	\$126.80
Scalpels, Sterile, Disposable, #11 Blade Individually wrapped. 10 blades/box.		6802A	1-4 bx	\$14.00
		6802A	5+ bx	\$11.00
Scalpel Handle, Stainless Steel	#4	5546B	each	\$12.60
	#8	5546	each	\$18.50

GROSSING & EMBEDDING

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
Scissors				
Standard & Unique Scissors				
Surgical Scissors: 17 cm Stainless Steel		5529	each	\$13.90
Scissors: 16.5 cm Curved Blunt/Blunt Tip Economy Grade		5530	each	\$4.30
Scissors: 16.5 cm Straight Blunt/Sharp Tip Economy Grade		5531	each	\$4.30
Mayo Scissors: 17 cm Straight Standard Grade		5532	each	\$25.40
Forman Scissors: 15 cm Angled	Standard Grade Stainless Steel	5533A	each	\$53.00
	Premium Grade Stainless Steel	5533B	each	\$76.00
	Tungsten Carbide	5533C	each	\$125.00
Cottle Scissors: 16 cm Angled	Standard Grade Stainless Steel	5534A	each	\$56.00
	Premium Grade Stainless Steel	5534B	each	\$77.00
	Tungsten Carbide	5534C	each	\$125.00
Flexible Probe with Eye, 13 cm, Stainless Steel Used to find and/or follow openings to vessels, ducts and various other structures.		5528	12/pk	\$39.60
Lighted Scissors				
Lighted Scissors Kit includes: one instrument, power pack and battery charger.	17 cm, Mayo Scissors Blunt/Blunt Tip	5525	kit	\$250.00
	17 cm, Surgical Scissors Blunt/Blunt Tip	5526	kit	\$250.00
Lighted Stitch Scissors Kit includes: one instrument, power pack and battery charger.	11.5 cm, Stitch Scissors Blunt/Blunt Tip	5527	kit	\$250.00
Needle Holder				
Needle Holder & Scissors: Olsen-Hager, 14 cm, Standard Grade		6310	each	\$16.50



PRODUCT NAME

SPECIFICATIONS

PART

VOLUME

PRICE

Rulers for Pathology

Opaque Metric Ruler

Flush cut at zero. Identical printing on both top and bottom edges. Printed identically on both sides.

15 cm

6651

10/pk

\$20.00

30 cm

6652

10/pk

\$25.00

Clear Metric Ruler

Clear, flexible plastic with metric on top and bottom.

30 cm

6650

10/pk

\$30.00

Tissue Dyes & Accessories

Squeezable Tissue Marking Dyes

Squeezable Dye Kit

18 squeezable dispenser. 3 each of black, blue, green, red, yellow and orange in foam reusable base. 4 ml dye/dispenser.

6900

\$60.00

Squeezable Dye Refills

15 pre-filled dispensers in convenient carton. 4 ml dye/dispenser.

Black

6901A

15/bx

\$50.00

Blue

6901B

15/bx

\$50.00

Green

6901C

15/bx

\$50.00

Red

6901D

15/bx

\$50.00

Yellow

6901E

15/bx

\$50.00

Orange

6901F

15/bx

\$50.00

Violet

6901G

15/bx

\$50.00

Tissue Marking Dye Bottles

Tissue Marking Dyes, 2 oz. Bottle

Black

6902A

each

\$32.00

Blue

6902B

each

\$32.00

Green

6902C

each

\$32.00

Red

6902D

each

\$32.00

Yellow

6902E

each

\$32.00

Orange

6902F

each

\$32.00

Violet

6902G

each

\$32.00

Tissue Marking Dyes, 8 oz. Bottle

Black

6903A

each

\$84.00

Blue

6903B

each

\$84.00

Green

6903C

each

\$84.00

Red

6903D

each

\$84.00

Yellow

6903E

each

\$84.00

Orange

6903F

each

\$84.00

Violet

6903G

each

\$84.00

Tissue Dyes with Racks

Tissue Marking Dyes, 2 oz., 7 color Kit w/Plastic Rack

7 - 2 oz. bottles. (Colors: black, blue, green, red, yellow, orange, violet)

6904P

each

\$225.00

GROSSING & EMBEDDING

PRODUCT NAME	SPECIFICATIONS	PART	VOLUME	PRICE
IPC Blue Biopsy Dye				
IPC Blue Biopsy Marking Dye Can be added directly to tissue processor. Does not interfere with H&E, special stains or IHC.	1 oz. Dropper Bottle	6899A	each	\$10.00
	2 oz. Flip-Top Bottle	6899B	each	\$17.50
	4 oz. Flip-Top Bottle	6899C	each	\$30.00
	16 oz. Flip-Top Bottle	6899E	each	\$85.00
Ink Aid				
Ink Aid Solution for setting tissue marking dye in special situations.	8 oz. Spray Bottle	6910	5/cs	\$52.50
Applicator Sticks				
Applicator Stick, Cotton Tip, 1/8" round		6905B	1,000/bx	\$15.50
Applicator Stick, Plain, 1/8" round		6905A	1,000/pk	\$11.00
Paraffin & Repellent				
PureAffin™ e Paraffin Melting Temp: 56-57°C		6532	8 x 1 kg bags/cs	\$95.00
PureAffin™ i Paraffin Melting Temp: 56-57°C With DMSO for enhanced filtration.		6533	8 x 1 kg bags/cs	\$95.00
PureAffin™ x Paraffin Melting Temp: 50-54°C		6534	8 x 1 kg bags/cs	\$95.00
Paraffin Repellent, 4 oz. An easy way to clean up wax from instruments and countertops. Also can be used as an effective mold release.		6535	each	\$35.00
Telescopic Glasses				
Telescopic Glasses 2X telescopes attached to glasses. Excellent for small detailed work.				
Plano (Clear) Lens - No Prescription		5631	each	\$399.00
Single Vision Lens - Customized to Your Prescription		5632S	each	\$449.00
Bifocal Lens - Customized to Your Prescription		5632B	each	\$449.00
Distance & Bifocal Lens - Customized to Your Prescription		5633	each	\$474.00
Mini High Intensity Light Assembly (attach to glasses)		5634	each	\$199.00



PRODUCT NAME

SPECIFICATIONS

PART

VOLUME

PRICE

Dissecting Needles

Dissecting Needles, Plastic

(Uses include: Separating and teasing tissue sections in water bath)

5220PL

12/pk

\$15.00

Hammer

Hammer, Premium Grade, Stainless Steel

26 cm, 250 gm

5548

each

\$135.00

Instrument Holder

Instrument Holder, Premium Grade, Stainless Steel

Ideal for holding knives, scalpels and other tools in the grossing area.

5808

each

\$179.00

PRODUCT NAME

PART

VOLUME

PRICE

Disinfectants

CaviCide & CaviCide1

CaviCide

Ready to use. Kills HBV, TB, HIV, VRE, MRSA and more.

24 oz. Spray bottle	ACP131024	each	\$13.90
	ACP131024	12/case	\$129.00
Gallon container	ACP131000	each	\$35.60
	ACP131000	4 x 1 gal	\$125.00
2.5 Gallon container	ACP131025	each	\$75.00
	ACP131025	2 x 2.5 gal	\$129.00

CaviCide1

Ready to use. Effective in killing 99.9% of bacteria, viruses and fungi in 1 minute.

24 oz. Spray bottle	ACP135024	12/case	\$144.00
Gallon container	ACP135000	4 x 1 gal	\$135.00
2.5 Gallon container	ACP135025	2 x 2.5 gal	\$157.00

CaviWipes & CaviWipes1

CaviWipes

Ready to use. Kills HBV, TB, HIV, VRE, MRSA and more.

160 Wipes/Canister	ACP131100	canister	\$14.00
	ACP131100CS	12 canisters/case	\$126.00

CaviWipes1

Ready to use. Effective in killing 99.9% of bacteria, viruses and fungi in 1 minute.

160 Wipes/Canister	ACP135100	12 canisters/case	\$140.00
--------------------	-----------	-------------------	----------

CaviWipes Canister Wall Bracket

Wall bracket to hold CaviWipes canister.

ACP131175	each	\$26.10
-----------	------	---------

Bleach Sprayer

Activate, 10:1 Bleach Sprayer

Dual-liquid sprayer automatically combines water and bleach.

Bleach sprayer, with refillable 11 oz. water cartridge. (Bleach cartridges not included.)	5101	each	\$30.50
Bleach cartridge refills, 11 oz. disposable	5102	12/case	\$53.20

Chlorine Test Strips (50/bottle)

5100	bottle	\$37.00
	6 bottles	\$204.00



PRODUCT NAME

PART

VOLUME

PRICE

Skin Usage

Vionex No Rinse Gel Antiseptic Handwash

For use when soap and water are not available.

18 oz. Pump bottle

ACP101628

each

\$15.00

Formalin Absorbing Pads

FAN Pad-GL

Used for cleaning up small spills, grossing on and using as a shelf liner to minimize formalin fumes.

Mini (11" x 10")

ABFAN-55

6 rolls/case

\$138.00

Small (6" x 14")

ABFAN-614

48/case

\$146.00

Large (9" x 21")

ABFAN-825

30/case

\$146.00

FAN Pad-GL Elite

Gross and dissect formalized tissue on. Non-woven blue laminate surface.

8" x 11"

ABFAN-811E

100/case

\$140.00

11" x 15"

ABFAN-1115E

100/case

\$230.00

FAN Pad Ultra

Gross and dissect formalized tissue on. Fabric surface.

8" x 11"

ABFAN-811

25 pads/pack

\$31.30

ABFAN-811

100 pads/case

\$119.00

11" x 15"

ABFAN-1115

25 pads/pack

\$51.20

ABFAN-1115

100 pads/case

\$199.00

Formalin Cube Containment

SafeCube Dispenser

Holds a 5 gallon formalin cube for dispensing and containing spills and leaks.

ABSCD-100A

each

\$419.00

Formalin Disposal

Formalex Green Formalin Neutralizer

Neutralizes formalin waste to a non-hazardous, non-toxic solution for sink disposal.

(Neutralization takes 8-12 hrs.)

ABFX-01

1 gal cube

\$35.00

ABFX-250

2.5 gal cube

\$70.00

ABFX-05

5 gal cube

\$125.00

ABFX-30

30 gal drum

\$648.00

ABFX-55

55 gal drum

\$998.00

PRODUCT NAME	PART	VOLUME	PRICE
Formalin Disposal			
Form-Zero Formalin Neutralizer			
Neutralizes formalin waste to a non-hazardous, non-toxic solution for sink disposal.			
(Neutralization takes 20-30 min.)	ABFZ550-12	12/case	\$87.00
	ABFZ550-32	32/case	\$222.40
Formalex Green Empty Spray Bottle	ABFX-32E	each	\$12.00
FRC-3™ Collection & Treatment Kit for Formalex Green	ABFRC-03K	kit	\$159.00
Includes: waste collection containers, graduated dispenser, test kit, tissue funnel, laminated instruction sheet, log sheet.			
Aldehyde Test Kit (100 test packets per kit)			
(Replacement Kit)	ABFRC-01K	kit	\$105.00
Replacement Aldehyde Waste Collection Container	ABFRC-250E	2.5 gal	\$12.00
Replacement Tissue Funnel with Screen	ABFX-FUN	each	\$12.00

Spill Kits & Components

Spill Response Kits			
All kits include: Safety goggles, gloves, waste bags and broom/pan.			
Polyform-F™ (For formaldehyde and glutaraldehyde)	ABSRK-PF3202	2 x 32 oz	\$69.00
	ABSRK-PF0102	2 x 1 gal	\$129.00
General Purpose (One of each: Acid, base, solvent, chlorine)	ABSRK-GP3204	4 x 32 oz	\$125.00
Custom (Mix & match your own custom combination. Limited to one Solvent Handler per kit. Add \$10.00 for each additional Solvent Handler/kit.)	ABSRK-GP-3204CO	4 x 32 oz	\$125.00
Spill Response Kits - Metal Cabinet			
All kits include: Safety goggles, gloves, waste bags and broom/pan.			
Custom (Mix and match your own custom combination. Limited to one Solvent Handler per kit. Add \$10.00 for each additional Solvent Handler/kit.)	ABSRKC-3205	5 x 32 oz	\$275.00
Acid Handler Spill Powder			
	ABAH-32	32 oz bottle	\$19.00
	ABAH-04	4 x 1 gal	\$179.00
Base Control Spill Powder			
	ABBC-32	32 oz bottle	\$19.00
	ABBC-04	4 x 1 gal	\$179.00
Chlorine Control Spill Powder			
	ABCC-32	32 oz bottle	\$25.00
	ABCC-04	4 x 1 gal	\$265.00
Solvent Handler Spill Powder			
	ABSH-32	32 oz bottle	\$39.90
	ABSH-04	4 x 1 gal	\$435.00



PRODUCT NAME	PART	VOLUME	PRICE
Polyform-F™ For formaldehyde and glutaraldehyde spills.	ABPF-32	32 oz each	\$19.00
	ABPF-04	4 x 1 gal/case	\$164.00
	ABPF-05	2 x 2.5 gal/case	\$164.00

Safety Glasses

Standard Safety Glasses

Standard Safety Glasses, Clear Lenses

Black	5635BLK	each	\$8.00
Pink	5635PNK	each	\$8.00
Blue	5635BLU	each	\$8.00

Bifocal Safety Glasses

Bifocal Safety Glasses, Clear Lenses

+1.0	5636A	each	\$17.00
+1.5	5636B	each	\$17.00
+2.0	5636C	each	\$17.00
+2.5	5636D	each	\$17.00
+3.0	5636E	each	\$17.00

Over Spectacle Safety Glasses

Over Spectacle Safety Glasses, Clear Lenses

Fits over regular prescription glasses. Small side windows and soft rubber temples. Ballistic Vo Rated.

5637	each	\$12.00
------	------	---------

Gloves Nitrile Powder Free

SafeSkin Nitrile Powder-Free Gloves

100 gloves/box

X-Small	UM026080	10 boxes/case	\$158.00
Small	UM026081	10 boxes/case	\$158.00
Medium	UM026082	10 boxes/case	\$158.00
Large	UM026083	10 boxes/case	\$158.00
X-Large	UM026084	10 boxes/case	\$158.00

PosiShield Nitrile Powder-Free Gloves

50 gloves/box

Small	5354S	box	\$13.00
Large	5354L	box	\$13.00

PRODUCT NAME

PART

VOLUME

PRICE

Air Purifier

GC Series

Treatment area is up to 800 sq. ft.

IQ Air Multi-Gas Filtration System, GC (Filters are included) Removes xylene, alcohol and formalin.	ST1CCUB0HGB	each	\$1,199.00
IQ Air Multi-Gas Cartridge, GC (Replacement filters)	ST102401600	4/pack	\$295.00
IQ Air VOC Filtration System, GC (Filters are included) Removes xylene, organic solvents.	ST1CDUB0HGB	each	\$1,199.00
IQ Air VOC Cartridge, GC (Replacement filters)	ST102401000	4/pack	\$295.00
IQ Air Chemisorber Filtration System, GC (Filters are included) Removes formalin.	ST1CBUB0HGB	each	\$1,249.00
IQ Air Chemisorber Cartridge, GC (Replacement filters)	ST102401800	4/pack	\$375.00
Pre-Filter HEPA, GC	ST102101200	each	\$78.00
Post Filter Sleeves, GC	ST102501000	4/pack	\$95.00

GCX Series

Treatment area is up to 1,200 sq. ft.

IQ Air Multi-Gas Filtration System, GCX (Filters are included) Removes xylene, alcohol and formalin.	ST1FCUB0LGB	each	\$2,199.00
IQ Air Multi-Gas Cartridge, GCX (Replacement filters)	ST102421600	4/pack	\$495.00
IQ Air VOC Filtration System, GCX (Filters are included) Removes xylene, organic solvents.	ST1FDUB0LGB	each	\$2,199.00
IQ Air VOC Cartridge, GCX (Replacement filters)	ST102421000	4/pack	\$495.00
IQ Air Chemisorber Filtration System, GCX (Filters are included) Removes formalin.	ST1FBUB0LGB	each	\$2,299.00
IQ Air Chemisorber Cartridge, GCX (Replacement filters)	ST102421800	4/pack	\$595.00
Pre-Filter HEPA, GCX	ST102141100	each	\$139.00
Post Filter Sleeves, GCX	ST102505000	4/pack	\$169.00

Sharps Containers

Sharps Container, Sharp Chute, Red 1.5 Quart
Autoclave and incinerator safe.

6813A	each	\$6.20
6813A	case of 18	\$107.10



PRODUCT INFORMATION

We are a distributor of Innovex Biosciences and Golden Bridge International Labs IHC product lines; Detection Kits, Chromogens, Antibodies and Specialty Products. Their part numbers and pricing are available on our website: www.newcomersupply.com

PRODUCT NAME	PART	VOLUME	PRICE
Buffers			
Phosphate Buffer 0.1M, pH 7.2-7.4	13313A	1 L	\$20.10
	13313B	1 gal	\$57.20
	13313C	10 L cube	\$78.00
	13313D	20 L cube	\$152.50
Phosphate Buffered Saline 0.01M, pH 7.2 (Use: IHC diluent, rinse and wash) Shelf Life: 6 mos.	133102A	1 L	\$14.00
	133102B	1 gal	\$18.80
	133102C	10 L cube	\$44.90
	133102D	20 L cube	\$74.90
Phosphate Buffered Saline 0.01M, pH 7.4 (Use: IHC diluent, rinse and wash) Shelf Life: 6 mos.	133132A	1 L	\$13.20
	133132B	1 gal	\$19.10
	133132C	10 L cube	\$42.10
	133132D	20 L cube	\$74.90
Phosphate Buffered Saline 0.05M, pH 7.4 (Use: Rinsing unfixed tissues and diluting antibodies, e.g., in IHC) Shelf Life: 6 mos.	133134A	1 L	\$13.20
	133134B	1 gal	\$19.90
	133134C	10 L cube	\$44.90
	133134D	20 L cube	\$72.20
Phosphate Buffered Saline 0.1M, pH 7.4 (Use: Rinsing unfixed tissues and diluting antibodies, e.g., in IHC) Shelf Life: 6 mos.	133104A	1 L	\$14.00
	133104B	1 gal	\$21.10
	133104C	10 L cube	\$39.60
	133104D	20 L cube	\$54.20
Tris Buffer 0.05M, pH 7.6 (Use: Working solution for IHC diluent, rinse and wash)	140301A	500 ml	\$24.90
	140301B	1 L	\$45.20
	140301C	1 gal	\$93.90
	140301D	10 L cube	Inquire
	140301E	20 L cube	Inquire
Tris Buffered Saline 0.05M, pH 7.6, 10X (Use: IHC wash solution; dilute 1:10 before use) Shelf Life 6 mos. Store at 4°C.	140304A	500 ml	\$18.30
	140304B	1 L	\$25.50
	140304C	1 gal	\$54.90
	140304D	10 L cube	\$167.10
	140304E	20 L cube	\$210.60



PRODUCT NAME

PART

VOLUME

PRICE

Retrieval Solutions**Citrate Buffer 0.01M, pH 6.0**

10355A	500 ml	\$18.10
10355B	1 L	\$24.90
10355C	1 gal	\$81.00

EDTA Buffer 0.001M, pH 8.0

1056A	500 ml	\$20.10
1056B	1 L	\$31.00
1056C	1 gal	\$73.20



2505 Parview Road • Middleton, WI 53562-2579 • 800-383-7799

www.newcomersupply.com

August 2015