

Safe Cube Dispenser™

The Ultimate workstation when filling containers with secondary containment.


Overview of Products:

Safe Cube Dispenser:

(Perfect for safe liquid storage and container filling)



20" Deep X 18.5" Wide X 17" High




Fan Pads:

(Super absorbent pads, specifically treated for absorbing & neutralizing formaldehyde, 10% formalin and glutaraldehyde solutions. Eliminates odors, too!)



Features of the Safe Cube Dispenser:

- 1) Made of durable inert plastic able to **contain a spill of up to 5 gallons** in main cabinet.






Distributed By:


NEWCOMER SUPPLY™
Your Histology Supplier

2505 Parview Road, Middleton, WI 53562 (800) 383 • 7799 (608) 831 • 7888 FAX (608) 831 • 0866 • www.newcomersupply.com


- 2) **Built in alarm** in both main cabinet and in dispenser tray to alert personnel immediately of spill.



- 3) Counter weights built into the Safe Cube main cabinet allow for container weight up to 40 lbs. in filling tray with **no risk of tipping**.



- 4) Internal sliding shelf makes it **effortless to position cube** for pouring and storage.



- 5) **Simple and easy to use!**

Ideal for Use In:

- Operating Room and Labor & Delivery
- Ambulatory Surgery
- Histology and Pathology Laboratories
- Clinical and Research Laboratories
- College and University Laboratories
- Morgue and Autopsy Suites

Distributed By:

NEWCOMER SUPPLY™
Your Histology Supplier

2505 Parview Road, Middleton, WI 53562 (800) 383 • 7799 (608) 831 • 7888 FAX (608) 831 • 0866 • www.newcomersupply.com

Ordering Information:

<u>Description</u>	<u>Volume</u>	<u>Order #</u>	<u>Price</u>
Safe Cube Dispenser	Each	ABSCD-100A	\$395.00
Fan Pad - GL Small	48/case	ABFAN-614	\$146.00
Fan Pad - GL Large	30/case	ABFAN-825	\$146.00

Please call Newcomer Supply if you have any further questions:

800-383-7799



Cubitainer® Combination Package and Cube® Insert are Registered Trademarks of the Hedwin Corporation, Baltimore, Maryland.

Distributed By:

NEWCOMER SUPPLY™
Your Histology Supplier

2505 Parview Road, Middleton, WI 53562 (800) 383 • 7799 (608) 831 • 7888 FAX (608) 831 • 0866 • www.newcomersupply.com