

# WARM UP YOUR SLIDES WITHOUT BURNING UP YOUR BUDGET!

It's simple, you need equipment that you can rely on every day of the week that won't break the bank. And those pesky little slides aren't gonna dry on their own very quickly! So let us help you out and save a little money in the process...

Both of these models are ideal for use in the lab for drying paraffin tissue sections on slides quickly, safely and with ease!

- **Easy to adjust and use**
- **Thermal heater ensures even heat transfer**
- **LED Temperature Display**
- **Temp. settings from room temp. to 70°C**



Part #6258

- 10" x 7" Surface
- Holds approx. 23 slides
- Approx. 10 lbs.
- Great for labs with limited bench space!

**\$395.00**



Part #6259

- 25" x 8" Surface
- Holds approx. 66 slides
- Approx. 13 lbs.
- Ideal for high volume labs!

**\$495.00**

The Pricing on these are so hot,  
I can cook a brat on them!...  
it's a Wisconsin thing, you know.

**For more info call:  
800-383-7799**



**NEWCOMER SUPPLY**

2505 Parview Road, Middleton, WI 53562 (608) 831-7888, Fax (608) 831-0866