


# Some Like It Hot™

## PARAFFIN DIS-

It's hard to keep cool when the pricing is **SO HOT!**



### HERE'S THE SCOOP ON THIS BEAUTY:

- ⇒ Large Capacity (2.5 Gallon) Table Top Dispenser
- ⇒ Built in Temperature Controller
- ⇒ Stainless Steel Interior for Easy Cleaning
- ⇒ Heat Tolerant Metal Faucet
- ⇒ Temperature Range: Rm Temp. to 72°C
- ⇒ Industrial Metal Outer Cabinet

<u>Product Description</u>	<u>Order #</u>	<u>Volume</u>	<u>Price</u>
Paraffin Dispenser, 2,5 Gallon Capacity	6527	Each	<b>\$699.00</b>

**WOW!**

Distributed By:

**NEWCOMER SUPPLY®**  
Your Histology Supplier

2505 Parview Road, Middleton, WI 53562 (800) 383 • 7799 (608) 831 • 7888 FAX (608) 831 • 0866 • [www.newcomersupply.com](http://www.newcomersupply.com)

Please visit our new website!