

WHAT IF THE FORMALIN CUBE SPILLS OR LEAKS?

- ☠ It could shut down your lab or entire department for an extended period of time
- ☠ You'll probably meet the Hazmat Team on less than optimal circumstances
- ☠ Your lab and anything near it will SMELL!
- ☠ The exposure will dramatically increase the health risk to you and co-workers
- ☠ To process the paperwork that it will generate costs BIG BUCK\$

No one can prevent spills or leaks, but we **CAN** contain and alert you to make your cleanup, expense & time minimal!


THE SAFE CUBE DISPENSER

(Unique Liquid Storage and Dispenser for 5 Gallon Cubes)



- ✓ Main cabinet can contain a spill up to 5 gallons!
- ✓ Filling tray can contain a spill up to 1/2 liter
- ✓ Built-in Alarm instantly notifies personnel of any spills
- ✓ Filling tray can hold up to 40 lbs. with no risk of tipping
- ✓ SIMPLE TO USE!!!


When not in use, the filling tray folds up for safe and convenient storage of formalin cubes.



(20" Deep X 18.5" Wide X 17" High)



We Have Spill Kits, Too!!!



(Side Panel View of Alarm)

Alarm sensors are located in both the Main Cabinet and Filling Tray



(Main Cabinet Sensor)



(Filling Tray Sensor)

Description	Volume	Order #	Price
Safe Cube Dispenser	Each	ABSCD-100A	\$395.00

Cubitainer® Combination Package and Cube® Insert are Registered Trademarks of the Hedwin Corporation, Baltimore, MD

Distributed by:



2505 Parview Road
Middleton, WI 53562

Phone (800) 383-7799

Fax (608) 831-0866

www.newcomersupply.com