

**RETRIEVE EASILY AND GENTLY AT  
60 C FOR 30 MINUTES OR 75 C FOR 15 MINUTES  
WITH...**

# UNI-TRIEVE

## LOW TEMPERATURE RETRIEVAL SOLUTION


UNIVERSALLY APPLICABLE TO \*ALL\* ANTIBODIES AND TISSUES  
ALSO FOR IHC, IMMUNOFLUORESCENCE, DNA & RNA PROBES

- NO Boiling
- NO High pH Buffer
- NO Citrate
- Contains NO Enzymes

YOUR CHOICE OF HEAT SOURCE:



Water Bath



Dry Bath



Hot Plate

**ONE RETRIEVAL SOLUTION FOR ALL ANTIBODIES AND TISSUES!**

**GENTLE! TISSUE DOES \*NOT\*  
FALL OFF SLIDES OR GET DAMAGED!**

**BACKGROUND-FREE**

**AVAILABLE NOW!**



Ordering Info:		
Product Number:	QTY:	Price:
NB325	1 liter	\$275



Distributed by:

**NEWCOMER SUPPLY**<sup>TM</sup>

2505 Parview Road  
Middleton, WI 53562

Phone (800) 383-7799  
Fax (608) 831-0866  
[www.newcomersupply.com](http://www.newcomersupply.com)

A Product of

**INNOVEX  
BIOSCIENCES**