


Affordable Lighted Water Bath!

With this NEW lighted water bath,
 "I just can't wait to get on the microtome
 again...cuttin' tissue with my friends..."


Sung to a Willie Nelson song...




New!




Overall Size: 14" x 14" x 4.5" (H)
 Glass Tray Size: 10" x 6.5" x 2" (H)




WATER BASIN IS LIGHTED FOR EXCELLENT VISIBILITY!



TEMPERATURE CONTROL RANGE: ROOM TEMP TO 75°C



LED DISPLAY PANEL FOR BOTH SET & ACTUAL TEMPERATURE



REMOVABLE GLASS BASIN

<u>Description</u>	<u>Volume</u>	<u>Order #</u>	<u>Price</u>
Tissue Water Bath	each	6950	\$600.00

Distributed By:



Your Histology Supplier

Middleton, WI • (800) 383 • 7799 • www.newcomersupply.com