

ADHESIVE MICROSCOPE SLIDE SAMPLE FORM

Name:

Company:

Address:

City: State: ZIP:

Phone: ext:

Email:

- Samples ship to U.S. and Canada only
- Return the sample form by:
email: dprine@newcomersupply.com
fax: 608-831-0866

Newsilane Part #: 5070

Adhesion: Very Strong
Application: Difficult and routine IHC. Difficult tissues with long incubations.
Unit: 72 slides/box, 1,440 slides/case
Price: \$58.00/box, \$874.00/case

MAS Part #: 5077

Adhesion: Very Strong
Surface: Hydrophilic
Application: Excellent for the most challenging tissues and applications.
Unit: 100 slides/box, 1,000 slides/case
Price: \$280.00/case

TOMO Part #: 5081

Adhesion: Strong
Surface: Hydrophilic
Application: For automated IHC platforms such as Ventana/Roche Benchmark and Leica Bond.
Unit: 100 slides/box, 1,000 slides/case
Price: \$280.00/case

MAS-GP Part #: 5079

Adhesion: Strong
Surface: Hydrophilic
Application: For all IHC applications, HIER, pH6 / pH9 protocols, routine histology, and frozen sections.
Unit: 100 slides/box, 1,000 slides/case
Price: \$280.00/case

APS Part #: 5078

Adhesion: Strong
Surface: Hydrophobic, treated to provide a permanent positive charge.
Application: IHC and routine histology.
Unit: 100 slides/box, 1,000 slides/case
Price: \$180.00/case

TRUBOND 380 Part #: 5080

Adhesion: Very Strong
Surface: Hydrophilic
Application: For HIER processed when tissue type requires extra strong adhesion, epitope enhancement and DNA probe procedures.
Unit: 100 slides/box, 1,000 slides/case
Price: \$280.00/case