

# ADHESIVE MICROSCOPE SLIDE SAMPLE FORM

Name:

Company:

Address:

City:  State:  ZIP:

Phone:  ext:

Email:

- Samples ship to U.S. and Canada only
- Return the sample form by:  
email: [dprine@newcomersupply.com](mailto:dprine@newcomersupply.com)  
fax: 608-831-0866

## NewSilane Part #: 5070

**Adhesion:** Very Strong  
**Application:** Difficult and routine IHC. Difficult tissues with long incubations.  
**Unit:** 72 slides/box, 1,440 slides/case  
**Price:** \$59.70/box, \$898.00/case

## MAS Part #: 5077

**Adhesion:** Very Strong  
**Surface:** Hydrophilic  
**Application:** Excellent for the most challenging tissues and applications.  
**Unit:** 100 slides/box, 1,000 slides/case  
**Price:** \$295.00/case

## TOMO Part #: 5081

**Adhesion:** Strong  
**Surface:** Hydrophilic  
**Application:** For automated IHC platforms such as Ventana/Roche Benchmark and Leica Bond.  
**Unit:** 100 slides/box, 1,000 slides/case  
**Price:** \$280.00/case

## MAS-GP Part #: 5079

**Adhesion:** Strong  
**Surface:** Hydrophilic  
**Application:** For all IHC applications, HIER, pH6 / pH9 protocols, routine histology, and frozen sections.  
**Unit:** 100 slides/box, 1,000 slides/case  
**Price:** \$280.00/case

## APS Part #: 5078

**Adhesion:** Strong  
**Surface:** Hydrophobic, treated to provide a permanent positive charge.  
**Application:** IHC and routine histology.  
**Unit:** 100 slides/box, 1,000 slides/case  
**Price:** \$180.00/case

## TRUBOND 380 Part #: 5080

**Adhesion:** Very Strong  
**Surface:** Hydrophilic  
**Application:** For HIER processed when tissue type requires extra strong adhesion, epitope enhancement and DNA probe procedures.  
**Unit:** 100 slides/box, 1,000 slides/case  
**Price:** \$280.00/case