

# Which Hematoxylin is best for you?



	<u>1201 Harris Modified</u>	<u>12013 Harris</u>	<u>1180A &amp; C Gill I</u>	<u>1180D &amp; F Gill II</u>	<u>1180G &amp; I Gill III</u>	<u>1202 Mayer Modified</u>	<u>1409 Weigert Set</u>
Method	Regressive	Regressive	Progressive	Progressive	Progressive	Progressive	Progressive
Formulation	Similar to Gill III	Harris	Single Strength	Double Strength	Triple Strength	No Chloral Hydrate	Iron Hematoxylin
Acetic Acid	Yes	No	Yes	Yes	Yes	Yes	NA
Does Mucin Stain?	No	No	No	Moderate	Yes	No	No
Does Glass Stain?	No	No	No	Moderate	Yes	No	No
Recommended Use	Good for automated stainers and high volume	Manual or automated staining	Cytology	Frozen & paraffin sections and IHC	Frozen sections and IHC darker staining stands up to high volume	Nuclear stain for non acidic special stains	Nuclear counterstain for acidic special stains
Shelf Life	2 years	2 years	1 year	1 year	1 year	2 years	2 years

All hematoxylin solutions manufactured by Newcomer Supply are mercury free and do not require daily filtering