

LABORATORY PEN SAMPLE FORM

Name:

Company:

Address:

City: State: ZIP:

Phone: ext:

Email:

- Samples ship to U.S. and Canada only
- Return the sample form by:
email: dprine@newcomersupply.com
fax: 608-831-0866



HistoTec™ Pen

- Fine Tip - 0.9 mm
- Black ink, dries quickly
- Withstands stains, xylene, and chemistry



Diamond® PerfectMark™ Pen

- Ultra Fine Tip - 0.75 mm
- Black ink, dries quickly
- Smearproof, waterproof, odorless, and light-resistant



KP Marker Plus Pen

- Ultra Fine Tip - 0.75 mm
- Black ink
- Writes easily on slides and cassettes



Moist Mark Plus™ Pen

- Fine Tip - 0.75 mm
- Black ink
- Writes without smearing on slides and cassettes