


SAFECUBE DISPENSER



Features:

- Main cabinet can contain a spill up to 5 gallons
- Filling tray can contain a spill up to 1/2 liter
- Built-in alarm notifies personnel of any spills
- Filling tray can hold up to 40 lbs. with no risk of tipping
- Internal sliding shelf easily positions cube for pouring and storage
- Non-slip rubber mat on bottom



Dimensions






Inner Dimensions



Outer Dimensions

Built in Sensors



Distributed By:

NEWCOMER SUPPLY
Your Histology Supplier